



User's Guide

For 32 Karat™
Software Version 4.0

149921AA

June 30, 2000

Beckman Coulter, Inc.
4300 N. Harbor Blvd., Fullerton, CA 92835
© Copyright 2000 Beckman Coulter, Inc.
Printed in U.S.A.


Safety Information


This section provides safety information and instructions for the hardware and accessories of the system.

Safety Symbols

The symbols displayed below and on the instrument are reminders that all safety instructions should be read and understood before installation, operation, maintenance or repair to this instrument is attempted.

When symbols are displayed in this manual, pay particular attention to the safety information associated with the symbol.

 WARNING	This icon accompanies text and/or other international symbols dealing with hazards to personnel. When present, it indicates that a potential hazard to your personal safety exists if information stated within the “ WARNING ” paragraph is not adhered to or procedures are executed incorrectly.
---	--

 CAUTION	This icon accompanies text and/or other international symbols dealing with potential damage to equipment. When present, it indicates that there is a potential danger of equipment damage, software program failure or that a loss of data may occur if information stated within the “ CAUTION ” paragraph is not adhered to or procedures are executed incorrectly.
---	--



HIGH VOLTAGE

Paragraphs marked by this symbol indicate that a potential hazard to your personal safety exists from a high voltage source. In this document, the “**WARNING**” icon will accompany this symbol.

To reduce risk of electrical shock, all devices employ a three-wire electrical cable and plug to connect the equipment to earth ground.


Ensure that the wall outlet receptacle is properly wired and earth grounded.


DO NOT use a three-to-two wire plug adapter.


DO NOT use a two wire extension cord or a two wire multiple-outlet power strip.


Disconnect power to the system before performing maintenance.


DO NOT remove any panels; panels should be removed only by qualified service personnel.


 <p>WARNING</p>	<p>BIOHAZARD</p> <p>Paragraphs marked by this symbol indicate that a potential hazard to your personal safety exists from a biological source. In this document, the “BIOHAZARD” icon will accompany this symbol.</p> <p>Normal operation of the instrument can involve the use of solvents and reagents that are toxic, flammable or biologically harmful.</p> <ul style="list-style-type: none">• Observe all precautionary information printed on the original solution containers.• Operate the system in the appropriate environment.• Take all necessary precautions when using pathology or toxic materials to prevent the generation of aerosols.• Observe all applicable precautionary procedures when using flammable solvents in or near the instrument.• Wear appropriate laboratory attire, e.g., safety glasses, gloves, lab coat and breathing apparatus, when working with hazardous materials.• Dispose of all waste solutions in a proper manner.
---	---


 <p>WARNING</p>	<p>SHARP OBJECTS/MOVING PARTS</p> <p>Paragraphs marked by this symbol indicate that a potential hazard to your personal safety exists from unblunted corners or other appendages on the outside or inside of the equipment and to moving parts. In this document, the “WARNING” icon will accompany this symbol.</p> <p>To avoid injury due to moving parts, observe the following:</p> <ul style="list-style-type: none">• Keep loose clothing and hair away from moving parts.• Never attempt to physically restrict movement of moving parts.
---	--

 <p>WARNING</p>	<p>HOT SURFACES</p> <p>Paragraphs marked by this symbol indicate that a potential hazard to your personal safety exists from heated surfaces or other appendages on the outside or inside of the equipment. In this document, the “WARNING” icon will accompany this symbol.</p>
---	--

 <p>WARNING</p>	<p>LASER LIGHT</p> <p>Paragraphs marked by this symbol indicate that a potential hazard to your personal safety exists from a laser source. In this document, the “WARNING” icon will accompany this symbol.</p> <p>Since laser light is not accessible to the user during normal operation of the system, the overall classification of these lasers is “Class 1”, i.e., “lasers which are safe under reasonably foreseeable conditions of operation.”</p> <p>To prevent eye damage from potentially harmful laser light, observe all safety warnings.</p>
---	---

 The symbol consists of a circle containing a vertical line that ends in three horizontal bars of decreasing width, representing a ground connection.	<p>PROTECTIVE EARTH OR GROUND TERMINAL</p> <p>This symbol identifies the location of the protective earth or ground terminal lug on the equipment.</p>
--	---

 The symbol is a simple rectangular frame with a vertical line extending upwards from the top center, representing a switch in the OFF position.	<p>OFF Position of Principal Power Switch</p> <p>This symbol graphically represents the equipment main power switch when it is in the OFF position.</p>
---	--

 The symbol is a simple rectangular frame with a vertical line extending downwards from the bottom center, representing a switch in the ON position.	<p>ON Position of Principal Power Switch</p> <p>This symbol graphically represents the equipment main power switch when it is in the ON position.</p>
---	--


Informations sur la sécurité


Cette section contient des informations et instructions relatives aux composants et accessoires du système. Elle se compose des sous-sections suivantes :


Symboles de sécurité


Les symboles illustrés ci-dessous et ceux apposés à l'instrument vous rappellent qu'il est crucial de lire et comprendre toutes les consignes de sécurité avant d'installer, d'utiliser, d'entretenir ou de réparer cet appareil.


Lorsque des symboles apparaissent dans ce mode d'emploi, respectez scrupuleusement les consignes qui les accompagnent.


 <p>AVERTISSEMENT</p>	<p>Cette icône accompagne le texte et/ou les autres symboles internationaux destinés à protéger le personnel contre les risques de blessures. Lorsque cette icône est présente, votre personnel risque de s'exposer à des dangers en cas de non-respect des consignes du paragraphe «Avvertissement» ou d'exécution incorrecte des procédures.</p>
--	---


 <p>ATTENTION</p>	<p>Cette icône accompagne le texte et/ou autres symboles internationaux associés aux risques de dommages matériels. Lorsque ce symbole est présent, il indique que l'équipement risque d'être endommagé, qu'une défaillance logicielle ou une perte de données risque de se produire en cas de non-respect des consignes figurant dans le paragraphe «Attention» ou d'exécution incorrecte des procédures.</p>
--	---


 <p>AVERTISSEMENT</p>	<p>HAUTE TENSION</p> <p>Les paragraphes marqués de ce symbole indiquent qu'une source de haute tension risque de compromettre la sécurité de votre personnel. Dans ce document, l'icône «AVERTISSEMENT» est toujours accompagné de ce symbole.</p> <p>Tous les dispositifs sont munis d'un cordon d'alimentation électrique à trois fils qui permet de les mettre à la terre et de diminuer le danger d'électrocution.</p> <p>Assurez-vous que la prise murale correspondante est correctement câblée et mise à la terre.</p> <p>N'utilisez AUCUN adaptateur permettant de passer de trois à deux fils.</p> <p>N'utilisez AUCUNE rallonge à deux fils ou une multiprise à deux fils.</p> <p>Coupez l'arrivée de courant au système avant d'effectuer une maintenance.</p> <p>N'enlevez AUCUN panneau ; ces derniers ne doivent être enlevés que par des techniciens de service qualifiés.</p>
---	--


 <p>AVERTISSEMENT</p>	<p>DANGER BIOLOGIQUE</p> <p>Les paragraphes marqués de ce symbole indiquent que des substances biologiques risquent de compromettre la sécurité de votre personnel. Dans ce document, l'icône «DANGER BIOLOGIQUE» est toujours accompagné de ce symbole.</p> <p>Des solvants et des réactifs toxiques, inflammables et présentent un danger biologique peuvent être utilisés pendant le fonctionnement normal de l'instrument.</p> <ul style="list-style-type: none">• Observez toutes les informations de précaution imprimées sur le récipient d'origine des solutions.• Utilisez le système dans un environnement adapté à l'application.• Lorsque vous utilisez des substances pathologiques ou toxiques, prenez toutes les précautions nécessaires afin d'éviter la création d'aérosols.• Observez toutes les mesures de précaution applicables quand vous utilisez des solvants inflammables sur l'appareil ou à proximité de celui-ci.• Portez des vêtements de protection pour laboratoire, tels que des lunettes, des gants, une blouse de laboratoire et un respirateur quand vous travaillez avec des substances dangereuses.• La mise au rebut des déchets de solutions doit se faire de manière appropriée.
---	---


 <p>AVERTISSEMENT</p>	<p>OBJETS TRANCHANTS/PIÈCES EN MOUVEMENT</p> <p>Les paragraphes marqués de ce symbole indiquent que des angles vifs ou d'autres pièces saillantes à l'extérieur et l'intérieur de l'appareil, ainsi que des pièces en mouvement risquent de compromettre la sécurité de votre personnel. Dans ce document, l'icône «AVERTISSEMENT» accompagne toujours ce symbole.</p> <p>Pour éviter d'être blessé par les pièces en mouvement, observez les mesures suivantes :</p> <ul style="list-style-type: none">• Éloignez les vêtements lâches et les cheveux des pièces en mouvement.• N'empêchez jamais le déplacement des pièces en mouvement.
---	---

 <p>AVERTISSEMENT</p>	<p>SURFACE CHAUDE</p> <p>Les paragraphes marqués de ce symbole indiquent que des surfaces chaudes ou d'autres pièces saillantes à l'extérieur et l'intérieur de l'appareil risquent de compromettre la sécurité de votre personnel. Dans ce document, l'icône «AVERTISSEMENT» est toujours accompagné de ce symbole.</p>
---	--

 <p>AVERTISSEMENT</p>	<p>FAISCEAU LASER</p> <p>Les paragraphes marqués de ce symbole indiquent qu'une source laser risque de compromettre la sécurité de votre personnel. Dans ce document, l'icône «AVERTISSEMENT» est toujours accompagné de ce symbole.</p> <p>Comme l'utilisateur ne peut pas accéder au laser pendant le fonctionnement normal du système, leur classification globale de ces lasers est «classe 1» ; en d'autres termes, «ces lasers ne présentent pas de danger dans des conditions de fonctionnement normales».</p> <p>Pour protéger les utilisateurs contre les rayons laser dangereux, tenez compte de tous les avertissements de sécurité.</p>
---	---

	<p>Terre de protection ou borne de mise à Terre</p> <p>Ce symbole désigne l'emplacement de la terre ou de la borne de mise à terre de l'équipement.</p>
---	--

	<p>Position HORS de l'interrupteur de courant principal</p> <p>Ce symbole représente graphiquement l'interrupteur principal de l'équipement en position HORS.</p>
---	--

	<p>Position EN de l'interrupteur de courant principal</p> <p>Ce symbole représente graphiquement l'interrupteur principal de l'équipement en position EN.</p>
---	--


Sicherheitsinformationen


Dieser Abschnitt enthält Sicherheitsinformationen sowie Anleitungen für System-Hardware und -Zubehör.


Sicherheitssymbole


Die unten dargestellten und am Gerät angebrachten Symbole sollen daran erinnern, daß alle Sicherheitsanweisungen zu lesen und zu verinnerlichen sind, bevor die Installation, Bedienung, Wartung oder Reparatur dieses Gerätes in Angriff genommen wird.


Wo im Handbuch Symbole erscheinen, sind die mit dem jeweiligen Symbol verbundenen Sicherheitsinformationen besonders zu beachten.


 <p>WARNUNG</p>	<p>Dieses Symbol begleitet Texte und/oder sonstige internationale Symbole, die auf Personengefährdung hinweisen. Das Vorliegen dieses Symbols weist auf mögliche Personengefährdungen hin, die auftreten können, falls die Verfahren in dem mit “WARNUNG” betitelten Abschnitt nicht korrekt ausgeführt werden.</p>
---	--


 <p>VORSICHT</p>	<p>Dieses Symbol begleitet Texte und/oder sonstige internationale Symbole, die auf potentielle Geräteschäden hinweisen. Das Vorliegen dieses Symbols weist auf mögliche Gefahrenquellen für das Gerät, auf Software-Programmausfälle oder auf Datenverluste hin, die auftreten können, falls die Verfahren in dem mit “VORSICHT” betitelten Abschnitt nicht korrekt ausgeführt werden.</p>
--	---




 <p>WARNUNG</p>	<p>HOCHSPANNUNG</p> <p>Mit diesem Symbol versehene Abschnitte weisen auf eine von einer Hochspannungsquelle ausgehende potentielle Personen-gefährdung hin. In diesem Dokument wird dieses Symbol vom “WARNUNG”-Symbol begleitet.</p> <p>Um das Risiko von elektrischen Stromschlägen zu reduzieren, sind alle Geräte mit dreiadrigen Kabeln und Steckern ausgestattet und können nur an geerdete Steckdosen angeschlossen werden.</p> <p>Sicherstellen, daß der Wandanschluß korrekt verkabelt und geerdet ist.</p> <p>AUF KEINEN FALL Steckeradapter für zweiadrige Steckdosen verwenden.</p> <p>AUF KEINEN FALL ein zweiadriges Verlängerungskabel oder eine zweiadrig verkabelte Mehrfachsteckdose verwenden.</p> <p>Vor Wartungsarbeiten die Stromversorgung des Systems unterbrechen.</p> <p>KEINERLEI Verkleidungen entfernen; Verkleidungen dürfen nur von dafür qualifizierten Kundendiensttechnikern entfernt werden.</p>
---	--

 <p>WARNUNG</p>	<p>BIOGEFÄHRDUNG</p> <p>Mit diesem Symbol versehene Abschnitte weisen auf eine von einer biologischen Substanz ausgehende potentielle Personengefährdung hin. In diesem Dokument wird dieses Symbol vom “BIOGEFÄHRDUNG”-Symbol begleitet.</p> <p>Der normale Betrieb des Systems kann den Einsatz von toxischen, brennbaren oder biogefährlichen Lösemitteln und Reagenzien umfassen.</p> <ul style="list-style-type: none">• Sämtliche auf den Original-Lösungsbehältern aufgedruckten Sicherheitsangaben beachten.• Das System in entsprechender Umgebung betreiben.• Beim Umgang mit pathologischen oder toxischen Substanzen alle notwendigen Vorsichtsmaßnahmen zur Vermeidung von Aerosolbildung ergreifen.• Beim Umgang mit brennbaren Lösemitteln am Gerät bzw. in dessen unmittelbarer Umgebung alle diesbezüglichen Vorsichtsmaßnahmen beachten.• Bei Arbeit mit Gefahrstoffen geeignete Laborkleidung tragen, wie bspw. Schutzbrille, Handschuhe, Laborkittel und Atemgerät.• Sämtliche Lösungsabfälle sind sachgerecht zu entsorgen.
---	---

 <p>WARNUNG</p>	<p>SCHARFKANTIGE OBJEKTE/ BEWEGLICHE TEILE</p> <p>Mit diesem Symbol versehene Abschnitte weisen auf eine von scharfen Kanten oder sonstigen hervorstehenden Objekten außerhalb oder innerhalb des Geräts sowie von beweglichen Teilen ausgehende potentielle Personengefährdung hin. In diesem Dokument wird dieses Symbol vom “WARNUNG”-Symbol begleitet.</p> <p>Zur Vermeidung von Verletzungen durch bewegliche Teile ist folgendes zu beachten:</p> <ul style="list-style-type: none">• Locker sitzende Kleidung und Haare von beweglichen Teilen fernhalten.• Niemals versuchen, die Bewegung beweglicher Teile von Hand einzuschränken.
---	--

 <p>WARNUNG</p>	<p>HEISSE FLÄCHEN</p> <p>Mit diesem Symbol versehene Abschnitte weisen auf eine von heißen Flächen oder sonstigen hervorstechenden Objekten außerhalb oder innerhalb des Geräts ausgehende potentielle Personengefährdung hin. In diesem Dokument wird dieses Symbol vom “WARNUNG”-Symbol begleitet.</p>
---	--

 <p>WARNUNG</p>	<p>LASERSTRAHLUNG</p> <p>Mit diesem Symbol versehene Abschnitte weisen auf eine von einer Laserlichtquelle ausgehende potentielle Personengefährdung hin. In diesem Dokument wird dieses Symbol vom “WARNUNG”-Symbol begleitet.</p> <p>Da der Anwender bei normalem Systembetrieb keinen Zugang zur Laserstrahlung hat, fallen diese Laser unter die allgemeine Klassifizierung “Klasse 1”, d.h., “Laser, die unter weitgehend vorhersehbaren Betriebsbedingungen keine Gefährdung darstellen”.</p> <p>Um Augenschäden durch die potentiell schädliche Laserstrahlung zu vermeiden, sind sämtliche Sicherheitshinweise zu beachten.</p>
---	---

	<p>SCHUTZERDE ODER MASSEKLEMME</p> <p>Dieses Symbol weist die Schutz Erde bzw. die Masseklemme am Gerät aus.</p>
	<p>AUS-Position des Hauptnetzschalters</p> <p>Durch dieses Symbol wird graphisch dargestellt, wenn der Hauptschalter für die Stromzufuhr auf "AUS" geschaltet ist.</p>
	<p>EIN-Position des Hauptnetzschalters</p> <p>Durch dieses Symbol wird graphisch dargestellt, wenn der Hauptschalter für die Stromzufuhr auf "EIN" geschaltet ist.</p>


Informazioni sulla sicurezza


Questa sezione contiene informazioni sulla sicurezza e istruzioni attinenti ai componenti meccanici e agli accessori del sistema.


Simboli di sicurezza


I simboli riportati sotto e apposti sullo strumento segnalano che prima di installare, usare, effettuare la manutenzione o riparare questo strumento occorre leggere attentamente tutte le istruzioni sulla sicurezza.


Ogniqualvolta si incontrano questi simboli nel manuale, leggere attentamente le informazioni sulla sicurezza da essi segnalate.


 <p>AVVERTENZA</p>	<p>Questo simbolo compare insieme al testo o ad altri simboli internazionali relativi ai rischi per il personale. Indica che se non si osserva scrupolosamente quanto segnalato o se si svolgono erroneamente procedure, si corre il pericolo di infortunio.</p>
--	--


 <p>ATTENZIONE</p>	<p>Questo simbolo compare insieme al testo o ad altri simboli internazionali relativi al rischio di danni all'apparecchiatura. Indica che se non si osserva scrupolosamente quanto segnalato o se si svolgono erroneamente procedure, si possono causare danni all'apparecchiatura, cessazione del funzionamento del programma software o perdita di dati.</p>
--	--


 <p>AVVERTENZA</p>	<p>ALTA TENSIONE</p> <p>I paragrafi contrassegnati con questo simbolo indicano che esiste il pericolo di folgorazione da una sorgente di alta tensione. Questo simbolo è accompagnato dal simbolo di “AVVERTENZA”.</p> <p>Per ridurre il rischio di scosse elettriche, tutti i dispositivi sono dotati di spina e cavo a tre conduttori per il collegamento dell’apparecchio con l’impianto di messa a terra.</p> <p>Accertarsi che la presa di corrente sia cablata adeguatamente e collegata all’impianto di messa a terra.</p> <p>NON utilizzare un adattatore per collegare una spina a tre conduttori con una a due conduttori.</p> <p>NON utilizzare un cavo di prolunga a due conduttori o una presa multipla a due conduttori.</p> <p>Prima di eseguire la manutenzione del sistema, scollegarlo dall’alimentazione.</p> <p>NON rimuovere i pannelli; devono essere rimossi esclusivamente dal personale di assistenza tecnica.</p>
--	---


 <p>AVVERTENZA</p>	<p>PERICOLO BIOLOGICO</p> <p>I paragrafi contrassegnati con questo simbolo indicano che esiste un pericolo per la propria incolumità a causa di una sorgente biologica. Questo simbolo è accompagnato dal simbolo di “AVVERTENZA”.</p> <p>Il funzionamento normale del sistema può comportare l’uso di solventi e reagenti tossici, infiammabili o biologicamente pericolosi.</p> <ul style="list-style-type: none">• Prendere tutte le precauzioni indicate sui contenitori originali della soluzione.• Fare funzionare il sistema in un ambiente adeguato.• Quando si usano materiali tossici o che possano causare condizioni patologiche, prendere tutte le precauzioni necessarie per prevenire la formazione di aerosol.• Quando si usano solventi infiammabili all’interno o in prossimità dello strumento, seguire tutte le procedure pertinenti per la sicurezza.• Quando si opera con materiali pericolosi, indossare un camice da laboratorio e usare occhiali e guanti di sicurezza e un respiratore.• Eliminare tutte le soluzioni di scarto in modo adeguato.
--	---


 <p>AVVERTENZA</p>	<p>OGGETTI TAGLIENTI/PARTI MOBILI</p> <p>I paragrafi contrassegnati con questo simbolo indicano che esiste un pericolo per la propria incolumità a causa di spigoli vivi o altre parti accessorie all’esterno o all’interno dell’apparecchio. Questo simbolo è accompagnato dal simbolo di “AVVERTENZA”.</p> <p>Per prevenire infortuni che possono essere causati dalle parti mobili, prendere queste precauzioni:</p> <ul style="list-style-type: none">• Tenere i capelli e indumenti larghi lontani dall’area delle parti mobili.• Non tentare mai di limitare il movimento delle parti mobili.
--	---

 <p>AVVERTENZA</p>	<p>SUPERFICIE AD ALTA TEMPERATURA</p> <p>I paragrafi contrassegnati con questo simbolo indicano che esiste un pericolo per la propria incolumità a causa di superfici o altre parti accessorie all'esterno o all'interno dell'apparecchio, ad alta temperatura. Questo simbolo è accompagnato dal simbolo di “AVVERTENZA”.</p>
--	--

 <p>AVVERTENZA</p>	<p>LUCE LASER</p> <p>I paragrafi contrassegnati con questo simbolo indicano che esiste un pericolo per la propria incolumità a causa di una sorgente laser. Questo simbolo è accompagnato dal simbolo di “AVVERTENZA”.</p> <p>Durante il funzionamento normale del sistema, la luce laser non è accessibile né osservabile, pertanto la classificazione laser complessiva di questi laser è “Classe 1”, ossia “laser che non mettono a rischio l'incolumità delle persone in condizioni di funzionamento ragionevolmente prevedibili”.</p> <p>Per proteggere gli utenti dai rischi della luce laser, osservare tutte le avvertenze sulla sicurezza.</p>
--	---

	<p>Collegamento con la terra o terminale di protezione</p> <p>Questo simbolo indica la posizione del collegamento con la terra o dello spinotto di protezione sull'apparecchio.</p>
---	--

	<p>Posizione OFF (SPENTO) dell'interruttore principale di alimentazione</p> <p>Questo simbolo indica il punto dell'interruttore principale di alimentazione su cui occorre premere per spegnere l'apparecchio.</p>
---	---

	<p>Posizione ON (ACCESO) dell'interruttore principale di alimentazione</p> <p>Questo simbolo indica il punto dell'interruttore principale di alimentazione su cui occorre premere per accendere l'apparecchio.</p>
---	---


Información sobre seguridad


En esta sección se proporciona información sobre la seguridad e instrucciones para la utilización del equipo y accesorios del sistema.


Símbolos de seguridad

Los símbolos ilustrados a continuación y que se encuentran sobre el instrumento recuerdan al usuario que debe leer y comprender plenamente todas las instrucciones relativas a la seguridad antes de la instalación, operación, mantenimiento o reparación de este instrumento.

Cuando vea esos símbolos en este manual, preste especial atención a la información de seguridad correspondiente a esos símbolos.


 <p>ADVERTENCIA</p>	<p>Este icono acompaña al texto y los símbolos internacionales que identifican riesgos potenciales para el personal. Su presencia significa que hay un riesgo potencial a la seguridad personal del operador si no se aplica la información incluida en el párrafo de la “ADVERTENCIA” o si no se ejecutan correctamente los procedimientos.</p>
---	---


 <p>PRECAUCIÓN</p>	<p>Este icono acompaña a texto y los símbolos internacionales que identifican riesgos potenciales de daños para el equipo. Su presencia indica un riesgo potencial de daños al equipo, fallos de programa o pérdida de datos si no se aplica la información incluida en el párrafo de la “PRECAUCIÓN” o si no se ejecutan correctamente los procedimientos.</p>
--	--


 <p>ADVERTENCIA</p>	<p>VOLTAJE ALTO VOLTAJE</p> <p>Los párrafos marcados con este símbolo contienen advertencias de riesgos potenciales a la seguridad personal debido a una fuente de alto voltaje. En este documento, junto a este símbolo se imprimirá el icono de “ADVERTENCIA”.</p> <p>Para reducir el riesgo de descargas eléctricas, todos los dispositivos están dotados de cable eléctrico de tres alambres y del enchufe adecuado para conectar el equipo a tomacorrientes con descarga a tierra.</p> <p>Verifique que el tomacorrientes de pared tenga una conexión y descarga a tierra adecuadas.</p> <p>NO use un adaptador de tres a dos alambres.</p> <p>NO use un cable de extensión de dos alambres ni una banda de enchufe múltiple de dos alambres.</p> <p>Antes de efectuar tareas de mantenimiento interrumpa el paso de corriente al instrumento.</p> <p>NO quite ninguno de los paneles; éstos sólo deben ser quitados por personal de servicio calificado.</p>
---	---


<div data-bbox="305 235 431 348" data-label="Image"> </div> <p data-bbox="261 369 479 396">ADVERTENCIA</p>	<p data-bbox="607 218 928 254">RIESGO BIOLÓGICO</p> <p data-bbox="607 275 1409 417">Los párrafos marcados con este símbolo indican riesgos potenciales a la seguridad personal debido a una fuente de sustancias biológicas. En este documento, junto a este símbolo se imprimirá el icono de “ADVERTENCIA”.</p> <p data-bbox="607 436 1409 541">La operación normal del sistema puede incluir uso de solventes y reactivos tóxicos, inflamables o que presentan un riesgo biológico.</p> <ul data-bbox="607 560 1409 1087" style="list-style-type: none"> • Observe todas las precauciones impresas en los recipientes originales de las soluciones. • Haga funcionar el sistema en el entorno adecuado. • Al utilizar materiales tóxicos o patológicos, tome todas las precauciones necesarias para impedir la generación de aerosoles. • Al utilizar solventes inflamables en el instrumento o cerca del instrumento, observe todas las precauciones aplicables. • Al trabajar con materiales peligrosos, use vestimenta de laboratorio adecuada, por ejemplo, gafas de protección, guantes, bata de laboratorio y dispositivo de respiración. • Procese todas las soluciones de desecho de la manera adecuada.
---	--


<div data-bbox="305 1236 431 1350" data-label="Image"> </div> <p data-bbox="261 1388 479 1415">ADVERTENCIA</p>	<p data-bbox="607 1205 1179 1241">OBJETOS AGUDOS/PIEZAS MÓVILES</p> <p data-bbox="607 1262 1409 1440">Los párrafos marcados con este símbolo indican riesgos potenciales a la seguridad personal debido a la existencia de esquinas u otras salientes agudas en el exterior o interior de este equipo. En este documento, junto a este símbolo se imprimirá el icono de “ADVERTENCIA”.</p> <p data-bbox="607 1459 1409 1528">Para evitar heridas provocadas por piezas móviles, observe las siguientes medidas:</p> <ul data-bbox="607 1547 1409 1690" style="list-style-type: none"> • Mantenga alejados de las piezas móviles la vestimenta y el cabello sueltos. • Nunca intente impedir físicamente el movimiento de las piezas móviles.
---	--

 <p>ADVERTENCIA</p>	<p>SUPERFICIES CALIENTES</p> <p>Los párrafos marcados con este símbolo indican riesgos potenciales a la seguridad personal debido a superficies o salientes a alta temperatura en el exterior o interior de este equipo. En este documento, junto a este símbolo se imprimirá el icono de “ADVERTENCIA”.</p>
---	--

 <p>ADVERTENCIA</p>	<p>LUZ DE LÁSER</p> <p>Los párrafos marcados con este símbolo indican riesgos potenciales a la seguridad personal debido a una fuente de luz de láser. En este documento, junto a este símbolo se imprimirá el icono de “ADVERTENCIA”.</p> <p>Durante el funcionamiento normal del sistema, el usuario no tiene acceso a la luz del láser. Por ello, el láser del sistema tiene la clasificación general de “Clase 1”, es decir, “dispositivos láser que son seguros en condiciones de funcionamiento razonablemente predecibles”.</p> <p>Para proteger a los usuarios de daños a los ojos por la luz potencialmente nociva del láser, observe todas las advertencias de seguridad.</p>
---	---

	<p>Terminal CON protección DE descarga a tierra</p> <p>Este símbolo identifica la ubicación en el equipo del terminal con protección de descarga a tierra.</p>
---	---

	<p>Posición de “desconectado” (“OFF“) del interruptor principal de corriente</p> <p>Este símbolo representa en forma gráfica el interruptor principal de corriente cuando está en la posición desconectado (OFF).</p>
---	--

	<p>Posición de “conectado” (“ON“) del interruptor principal de corriente</p> <p>Este símbolo representa en forma gráfica el interruptor principal de corriente cuando está en la posición conectado (ON).</p>
---	--


安全性に関する情報


本セクションではシステムのハードウェア並びにアクセサリの安全性に関する情報及び指示を記載しています。

マーク


以下のマークが装置上に貼付されている場合、インストール、操作、メンテナンス、あるいは修理を行う前に安全に関する説明を読み、理解していることを再確認してください。


また、本書中にこれらのマークが使用されている場合は特に安全に関する説明を注意してお読みください。


 <p>要注意</p>	<p>このマークが文中や他の国際的に使用されているマークと共に使用されている場合、人体に害があることを示しています。「要注意」と記されている内容を守らなかったり誤って実施した場合、使用者に害を及ぼす可能性があることを意味します。</p>
--	--


 <p>注意</p>	<p>このマークが文中や国際的に広く使用されている他の国際マークと共に使用されている場合、装置に害があることを示しています。「注意」と記されている内容を守らなかったり誤って実施した場合、装置に損傷を与えたり、ソフトウェアにプログラム エラーが生じたり、またデータが消失する可能性があることを意味します。</p>
---	---


<div data-bbox="310 306 435 420" data-label="Image"> </div> <div data-bbox="324 436 417 470" data-label="Text"> <p>要注意</p> </div>	<div data-bbox="600 289 704 325" data-label="Section-Header"> <h3>高電圧</h3> </div> <div data-bbox="600 346 1417 462" data-label="Text"> <p>このマークは高電圧の危険性を意味し、使用者に害を及ぼす可能性があることを示しています。本書では「要注意」のアイコンと共に使用されています。</p> </div> <div data-bbox="600 531 1417 604" data-label="Text"> <p>感電防止のために、全ての装置は 3 線式の電気コードを用い接地された装置に接続してください。</p> </div> <div data-bbox="600 672 1417 745" data-label="Text"> <p>壁面のコンセントは正しく配線され接地されていることを確認してください。</p> </div> <div data-bbox="600 814 1287 850" data-label="Text"> <p>3 線から 2 線用のアダプタは使用しないでください。</p> </div> <div data-bbox="600 917 1417 991" data-label="Text"> <p>2 線用の延長線または 2 線用のテーブルタップは使用しないでください。</p> </div> <div data-bbox="600 1060 1385 1096" data-label="Text"> <p>メンテナンスを行う前にシステムの電源を切ってください。</p> </div> <div data-bbox="600 1163 1417 1236" data-label="Text"> <p>いかなるパネルも取り外さないでください。パネルは弊社サービス員が取り外します。</p> </div>
---	--


 <p>要注意</p>	<p>生物学上の危険性</p> <p>このマークは生物由来による人体への害の可能性があることを示しています。本書では「要注意」のアイコンと共に使用されています。</p> <p>通常の操作でも有毒、可燃性、あるいは生物学上有害な溶媒や試薬を使用することがあります。</p> <ul style="list-style-type: none">・ 容器に記載された注意事項をよくお読みください。・ 適切な環境下で装置を操作してください。・ 病理上または有毒物質を使用する際はエアゾルの発生を防ぐための予防措置をとってください。・ 装置のそばで可燃性の溶媒を使用する際は該当する安全手順に従ってください。・ 有害物質を取り扱う際は適切な保護眼鏡、手袋、上着等を着用してください。・ 廃液は適切な方法で廃棄してください。
--	---

 <p>要注意</p>	<p>鋭利物 / 稼動部分</p> <p>このマークは鋭利な角や装置内外の付属物、及び稼動部分による人体への害の可能性があることを示しています。本書では「要注意」のアイコンと共に使用されています。</p> <p>可動部分による傷害を防ぐために以下の点を確認してください。</p> <ul style="list-style-type: none">・ ゆったりした衣服や長髪は稼動部分から離してください。・ 稼動部分の動きを手で止めようとしないでください。
--	---

 <p>要注意</p>	<p>高温箇所</p> <p>このマークは高温箇所や装置内外の付属物による人体への害の可能性を示しています。本書では「要注意」のアイコンと共に使用されています。</p>
--	---

 <p>要注意</p>	<p>レーザー光</p> <p>このマークはレーザー光線による人体への害の可能性を示しています。本書では「要注意」のアイコンと共に使用されています。</p> <p>通常の操作ではレーザー光はユーザの目に触れることはありません。そのため、システム全体としてのレーザーの分類は「クラス 1」、すなわち「操作上無理なく予測できる状態の下での安全なレーザー」となっています。</p> <p>レーザー光の害から目を守るために、すべての注意書きを確認してください。</p>
--	---

	<p>接地用端末</p> <p>このマークは装置に接地用端末があることを示しています。</p>
---	---

	<p>電源スイッチ OFF 状態</p> <p>このマークは装置の主要電源スイッチが OFF の状態になっていることを図により示しています。</p>
---	--


	<p>電源スイッチ ON 状態</p> <p>このマークは装置の主要電源スイッチが ON の状態になっていることを図により示しています。</p>
---	--

Table of Contents

Instructions for using the Data Samples	vi
I. This folder contains the following files:	vi
II. Purpose of these files.....	vi
III. Using the files	vi
IV. IMPORTANT NOTICE: Duplicate these files before use	vii
Introduction to CE	
Overview	1-1
General Description	1-1
Common Modes of Electrophoresis	1-2
Capillary Zone Electrophoresis (CZE)	1-2
Micellar Electrokinetic Capillary Chromatography (MEKC)	1-2
Capillary Gel Electrophoresis (CGE)	1-2
Capillary Isoelectric Focusing (cIEF).....	1-2
Capillary Electrochromatography (CEC)	1-3
System Overview	
The Sample Handling System.....	2-2
The Capillary Cartridge	2-5
The Syringe Pump	2-6
The High Voltage Power Supply.....	2-6
LED Indicators.....	2-7
Cartridge and Tray Cover Interlocks.....	2-8
Opening the tray cover.....	2-8
Opening the cartridge cover.....	2-8
UV Detector Optics	2-9
The Photo Diode Array (PDA) Detector.....	2-11
Laser Induced Fluorescence Detector.....	2-13
The Laser Modules	2-15
488 nm Laser Module.....	2-16
635 nm Laser Module.....	2-17
Using Other Lasers with the LIF Detector.....	2-17
Starting 32 Karat Software	
Tool Bar	3-2
Menu Bar	3-2

Creating and Configuring an Instrument	3-3
Creating a New Instrument	3-3
Configuring the New Instrument	3-4
Auto Configuration	3-7
Configuration Options	3-9
Analysis Options.....	3-10
Instrument Options	3-11
Manual Configuration.....	3-12
Manual Configuration of the UV Detector	3-12
Configuration Options	3-13
Manual Configuration of the PDA and LIF Detectors.....	3-13

Direct Control

Introduction.....	4-1
Direct Control Window.....	4-2
Exercise: Conditioning the Capillary.....	4-5

Method Editing

Creating a Method	5-1
Introduction.....	5-1
Initial Conditions Tab	5-1
Auxiliary data channels	5-2
Mobility channels	5-2
Temperature	5-2
Peak detect parameters.....	5-3
Trigger settings	5-3
Inlet and Outlet trays.....	5-3
Detector Initial Conditions.....	5-3
UV Detector Initial Conditions.....	5-4
Electropherogram channel	5-4
Filter.....	5-5
Relay 1 and Relay 2.....	5-5
Absorbance Signal	5-5
PDA Detector Initial Conditions.....	5-6
Electropherogram Scan Data	5-6
Filter.....	5-6
Relay 1 and Relay 2.....	5-6
Reference Channel.....	5-6
Absorbance Signal	5-6
Electropherogram Channel Data.....	5-7

Multi-Egram	5-7
LIF Detector Initial Conditions.....	5-8
Electropherogram Channel 1 and 2	5-8
Signal	5-9
Data Rate.....	5-9
Relay 1 and 2	5-9
Time Program	5-10
Separate.....	5-12
Rinse	5-14
Inject	5-15
Autozero.....	5-16
Stop Data.....	5-16
End	5-16

Running the System

Single Run.....	6-1
Materials needed:	6-1
Programming a Sequence	6-5
Running the Sequence	6-11

Integration

Description.....	7-1
Integrating Data	7-2
Optimizing Integration.....	7-4
Integration Parameters	7-5
Width.....	7-5
Threshold	7-7
Integrating a Peak Cluster.....	7-8
Integration Off.....	7-8
Integration Results and Reports	7-9

Calibration

Introduction.....	8-1
Qualitative Analysis	8-1
Quantitative Analysis.....	8-1
Creating a Calibration.....	8-2
Developing a Calibration Method.....	8-2
Generating the Calibration Curve	8-4
Analyzing Unknowns.....	8-8

Qualitative Analysis

Data Analysis and Reporting9-1

Using Mobility

Mobility Markers10-2
 Assigning Mobility Values.....10-2
 Mobility Plot.....10-4
Demonstration: Using the Mobility Features.....10-5
Peak Identification by Mobility10-8

Creating Reports

Introduction.....11-1
Creating a Method Custom Report11-2
 Modifying a Report Template.....11-2
 Saving a Report as a Template.....11-4

Appendix 1

Additional Resources Appendix 1-1
 32 Karat Online Help Appendix 1-1
 Installation and Maintenance Manual Appendix 1-1
 Beckman Coulter Service Appendix 1-1
 Help by e-mail Appendix 1-1
 Other references Appendix 1-2

Appendix 2

CE Calculations Appendix 2-1
 Ohm's Law Appendix 2-1
 Corrected Migration Time Appendix 2-1
 Corrected Area..... Appendix 2-1
 Mobility Appendix 2-2

Appendix 3

Understanding PDA Data Appendix 3-1
 Mixed View w/ 3D option..... Appendix 3-2
 Contour plot Appendix 3-2
 Spectrum Appendix 3-3
 Electropherogram..... Appendix 3-3
 3D Plot Appendix 3-3

Appendix 4

System Administration Appendix 4-1

Preface



The 32 Karat™ User's Guide is a self-paced introduction to the operation of the Beckman Coulter, Inc. P/ACE MDQ and the 32 Karat Software version 4. It contains a brief introduction to the technique of capillary electrophoresis; an overview of the MDQ hardware; tutorials for the most frequently used software features; and an appendix of useful information.

The 32 Karat User's Guide has been designed to assist users who are new to CE or the 32 Karat Software in beginning to use both the software and the instrument. It is largely built on exercises that are designed to take you step-by-step to a point where you can begin to create your own methods, generate and analyze data, and report results. It does not attempt to cover every possible feature in the system or every aspect of capillary electrophoresis.

The on-instrument activities contain specific instructions for users with ultraviolet, photodiode array, and laser induced fluorescent detector systems. Follow the instructions for the detector(s) installed on your system.

The activities in this manual are progressive. Each builds on the previous ones. To get the most from this manual, work through all the activities, in order.

The activities that teach data reprocessing, calibration, and reporting, are all based on data collected with an ultraviolet detector system. The steps in these procedures are very similar regardless of the detector type used, so these exercises are important for users of all detector types. Data files for the post-run exercises are installed with the 32 Karat Software. If your instrument does not have a UV detector, you can manually configure a virtual UV instrument in the software so that you can do these exercises.

This manual was created by a team of technical writers, scientists, and engineers at Beckman Coulter, Inc. We value your feedback. Please let us know how we can improve future editions of this manual. Send your comments to: 32karat@beckmancoulter.com, and include User Manual in the subject line. Questions concerning the exercises can be sent to the same address.

Fullerton, California

June, 2000

Instructions for using the Data Samples

In some of the exercises in this manual you are asked to open specific data and method files. These are located in the Data Samples folder, which is a sub-directory to the 32Karat folder that was created when your software was installed. The following information is also contained in the "Readme.doc" file located in the Data Samples folder. It is IMPORTANT to read and understand the following before using the sample files.

I. This folder contains the following files:

- Data Sample 1.dat
- Data Sample 2.dat
- Data Sample 3.dat
- Level 1.dat
- Level 2.dat
- Level 3.dat
- Level 4.dat
- Level 5.dat
- PDA Data Sample.dat
- Calibrate.met
- Quantitate.met

II. Purpose of these files

The Data Samples folder contains data and method files to supplement the exercises in the User's Guide for the 32 Karat Software. The files were created solely for the purposes of this tutorial. They are intended as training tools only, and are not to be used for system validation or any other purposes.

III. Using the files

These files are referred to by name in the 32 Karat User's Guide. In performing some of these exercises, you will be asked to open one or more of these files. By using these files, the images you see on the computer screen and the images in the training manual will be the same.

IV. IMPORTANT NOTICE: Duplicate these files before use

Processing these data files and methods in the 32 Karat Software will cause changes to the files. To keep the original files unaltered and available for later use in these exercises, please follow the steps below.

1. Open Windows NT Explorer. To open this utility, right click on the Start button in the Windows NT toolbar and select Explore.
2. Navigate to the Data Samples folder. It is a sub-directory of the 32Karat folder.
3. In the left-hand pane of the Explorer window, click on the Data Samples folder to open it.
4. From the menu bar at the top of the Explorer window select File | New | Folder. A new folder will be added to the Data Samples folder. The “New Folder” will be highlighted, which indicates that the name may be edited. Type a new name. The name of the individual who will be using the tutorial exercises in the User’s Guide is a good choice.
5. Select all of the data files listed above, but do not select the new folder. There are several ways to do a multiple file select. Use any one of the following:
 - Select the first file in the list by clicking on it. Hold down the Shift key, and click on the last file in the list. All the files between these two files will be selected.
 - Select a file by clicking on it. Hold down the Ctrl key, and click each individual file in the list.
 - Using the mouse, you can click and drag a box around the files to be selected.
6. When all the files are selected (highlighted), select Edit | Copy from the menu bar (or use the keyboard shortcut Ctrl-C). Double click on the new folder to open it. Then select Edit | Paste from the menu bar (or use the keyboard shortcut Ctrl-V).
7. You have now created a copy of the sample files in a new folder. When the manual asks you to open a file, open the copy you have just created, and not the original file. **DO NOT USE THE ORIGINAL FILES!** Save all methods and sequences you create as part of the tutorial in this new folder.
8. A backup set of files is included in case the originals become altered. These are stored in the file “Backup Files.exe.” Running this executable file will restore the original set of sample files. When running the backup files program, you may be prompted to over-write existing files. Select “Yes.” The contents of the Data Samples file will be replaced by this action, including any changes you have made to the original set of data and method files.

Section 1-Introduction to CE

Overview

This section describes the common modes of capillary electrophoresis and how these processes are performed by the P/ACE™ MDQ instrument using the Ultraviolet (UV), Photodiode Array (PDA) and Laser Induced Fluorescence (LIF) detectors.

General Description

The P/ACE MDQ instrument separates sample components within a fused-silica capillary, using one of several modes of electrophoresis. All of these modes are generally referred to as Capillary Electrophoresis (CE).

In the P/ACE MDQ, sample is injected into the capillary using either pressure, vacuum, or voltage. Under the influence of high voltage, sample components migrate differentially through the capillary.

There are two basic methods for detecting samples using CE. The first is the absorbance of light. This method is utilized in UV and PDA detection. As these components pass a window in the capillary, a single wavelength UV detector or a multiwavelength photodiode array detector (PDA) measures absorbance and transmits the signal to the computer. The signal can also be transmitted to an external recorder, integrator, or data system through an analog output. The signal may be plotted graphically in the form of an electropherogram and analyzed.

The second method is to induce the samples to fluorescence and measure the emitted light. This is done by Laser Induced Fluorescence (LIF). Substances in the capillary which fluoresce at the laser wavelength are detected. The LIF Detector measures and records this fluorescence, which appears as a peak on the computer screen or printed electropherogram.

The P/ACE MDQ instrument may be used to separate many different kinds of samples, including peptides, proteins, nucleic acids, ions, enantiomers, and pharmaceuticals. CE has very low sample requirements relative to other analytical techniques (5 to 30 microliters, with the actual injection volume being typically between 5 and 50 nanoliters). It provides a complimentary alternative to other separations techniques such as chromatography.

Common Modes of Electrophoresis

Capillary Zone Electrophoresis (CZE)

When voltage is applied to an uncoated fused-silica capillary tube filled with a uniform electrolyte solution, separation of charged species occurs due to differential migration in the electrical field. This process, in which charged particles in solution migrate toward an electrode with opposite charge, is called electrophoresis. If the inner wall of the capillary is also charged, the fluid in the capillary will begin to flow toward the electrode that has the same charge as the capillary wall. This bulk movement of fluid is termed electroosmotic flow or EOF.

Because EOF is usually of higher magnitude than electrophoresis, analytes of both positive and negative charge will ultimately be carried in the same direction, although at different rates. Electrophoresis will separate like-charged but otherwise distinct analytes from one another. In this manner, both positive and negative analyte molecules can be detected as they are swept past the detector. The direction of EOF can be changed by reversing the charge on the electrodes, or by changing the charge on the wall of the capillary by chemical means.

The amount of time required for sample molecules to migrate to the detector depends on the length of the capillary, the electrophoretic mobilities of the particular sample molecules, the specific electrolyte used, the magnitude of the EOF, and the applied voltage. Molecules having different electrophoretic mobilities will be detected at different times. Other factors to be considered are the size, shape, and charges of the particles, the electrolyte concentration, the pH of the separation buffer, and the dimensions of the capillary.

Micellar Electrokinetic Capillary Chromatography (MEKC)

In free-zone capillary electrophoresis, neutral molecules travel as a single band. To separate these molecules, micellar additives such as sodium dodecyl sulfate (SDS) can be added to the electrolyte solution. Neutral molecules with differing affinities for the charged micelles will separate during electrophoresis under these conditions.

Capillary Gel Electrophoresis (CGE)

In the techniques described above, the analyte molecules are moving in a low viscosity liquid. By increasing the viscosity of the matrix to the point where the migration of molecules is physically impeded, it is possible to separate different species based on size and shape, as well as charge. A variety of gel solutions may be employed in this technique.

Capillary Isoelectric Focusing (cIEF)

In IEF, a mixture of special buffers called ampholytes is used to create a pH gradient along the length of the capillary. During separation, molecules will migrate to the point in the pH gradient at which they have no net charge. At the end of the process, the various analytes are located in discrete zones in the capillary. In a separate mobilization step, the bands are moved past the detector for analysis. This technique is most often used for protein and peptide analysis.


Capillary Electrochromatography (CEC)

CEC uses a capillary packed with porous particles similar to those used in HPLC. Driven by EOF, analytes partition between mobile and stationary phases as they migrate down the capillary. Under some conditions, electrophoresis may also participate in the separation. CEC is primarily employed to separate small molecules such as drugs.

Section 2-System Overview

The main components of the P/ACE MDQ instrument include trays that hold vials of sample, buffer, and other solutions, an interface block, a high-voltage power supply and electrodes, a source optics module and detector, temperature control hardware, and a sample injection mechanism.

Figure 1 P/ACE MDQ instrument



- | | |
|--|------------------------|
| 1 Source Optics Module with D ₂ Lamp inside | 9 Power Switch |
| 2 Insertion Bar | 10 Outlet Buffer Tray |
| 3 Interface Block | 11 Outlet Sample Tray |
| 4 High Voltage Electrode | 12 Grounded Electrode |
| 5 High Voltage Power Supply (inside) | 13 Fiber Optic Cable |
| 6 Inlet Sample Tray | 14 Detector |
| 7 Inlet Buffer Tray | 15 Capillary Cartridge |
| 8 Coolant Fill Ports | |


The main power switch is located on the lower right side of the front of the instrument. All connections for external system components are located on the upper left side panel of the instrument, except for AC inlet and the fuse holder. Three fans provide cooling airflow for internal system components; air is exhausted through vents at the side and rear of the instrument. Maintain proper clearances to ensure that the vents are not blocked (see the P/ACE MDQ Installation and Maintenance Manual).

The Sample Handling System

The sample handling system holds four trays; two sample trays (inlet and outlet), and two buffer trays. The sample trays are primarily used for samples; the buffer trays hold the other solutions required for electrophoresis (e.g. buffer and rinse solutions). The trays are arranged on two parallel tracks. Under normal operating conditions, the trays on the left are referred to as the inlet trays for sample and buffer; the trays on the right are referred to as the outlet trays for sample and buffer.


Each buffer tray holds either 36-2mL vials, or a large reservoir tray; the large reservoir tray has two 25 mL buffer reservoirs and holders for ten 2mL vials. Sample trays hold either 48 vials or a 96 well plate. Each tray is numbered from the front to the back, beginning with the number 1, and lettered from left to right, beginning with the letter A (see Figure 2).


Figure 2 Sample trays



901018L.AI

- 1 Buffer Tray
- 2 48 Vial Sample Tray
- 3 96 Position Sample Tray
- 4 Large Reservoir Tray


 <p>WARNING</p>	<p>The P/ACE MDQ is not intended to accommodate volatile materials in 96-well plates. Volatile solvents may release hazardous or flammable vapors leading to a risk of fire or explosion. Solvent vapors may also damage the instrument. Do not use volatile solvents in 96-well plates.</p>
--	--

 <p>WARNING</p>	<p>The 2 mL vials are pressurized during rinse and separation-with-pressure events. To reduce the risk of breakage and expelled glass particles, use only Beckman vials P/N 144980, and inspect each vial for damage prior to use. Do not use any vial that appears to be cracked or damaged in any way. Wear safety glasses when opening the sample cover while the vials are pressurized.</p>
--	---

The Capillary Cartridge

The CE capillary is installed in a cartridge. The cartridge design protects the capillary, provides a path for liquid coolant, simplifies installation into the instrument, and aligns the detection window in the optics. The components of the cartridge are shown in Figure 8.

Figure 3 The capillary cartridge



901011L.AL

1. Coolant tubing with capillary inside
2. Capillary

3. Detector Window and Aperture -
UV and PDA Detectors
4. Detector Window and Aperture - LIF Detector

The detection window is an area of the capillary where the polyimide coating has been removed to reveal the transparent fused-silica. This area of the capillary is placed in a part of the cartridge containing a plug that connects the window to the optical system. One type of plug is used for UV and PDA detectors, and a second type is used for LIF detectors. The procedures for installing a capillary in a cartridge are covered in the Installation and Maintenance Manual. In addition to the cartridge shown here, special cartridges are used with the optional External Detector Adapter (see the Installation and Maintenance Manual for information on these cartridges.)

The capillary temperature is controlled with an inert liquid that circulates through the cartridge. Temperature is controlled in a range from 10°C below ambient (with a minimum of 15°C) to 60°C. Coolant flows through the cartridge via two openings in the bottom of the housing (located between the ends of the capillary.) This fluid removes the heat generated by electrophoresis.

The Syringe Pump

P/ACE MDQ is capable of generating pressures using an internal pump mechanism. This pump can supply 0.1 to 25 psi to perform pressure injections or low pressure mobilizations. It can also apply up to 100 psi to move fluids through the capillary. Vacuum injections can be performed from 0.1 to 5.0 psi. Pressure can be applied to both ends of the capillary simultaneously to facilitate methods such as capillary electrochromatography.


The High Voltage Power Supply

The high voltage power supply is capable of delivering a maximum of 30 kV with a maximum current of 300 μA . The voltage range is from 1 to 30 kV, selectable in 100 V increments. The polarity is selectable from the software. The current is selectable from 3.0 to 300 μA in 0.1 μA increments. The software allows the user to select current, voltage or power operation. During operation, the system will ramp the voltage or current up to the programmed value. Limits for voltage, current and power may also be entered to protect the capillary. For example, if the user programs a voltage setting for 30 kV but the setting for current is only 3.0 μA , the system reaches the limit set for current before reaching the voltage setting, and regulates voltage to maintain that current.

LED Indicators

The front panel of the instrument contains LED indicators for Power, UV, and High Voltage.

Figure 4 LED Indicators



1. High Voltage Indicators
2. UV Lamp On Indicator

3. Instrument Power Indicator

Cartridge and Tray Cover Interlocks

The double-hinged doors of the P/ACE MDQ have been designed with interlock sensors that will prevent unsafe access to the interior of the instrument. The first (lower) door is called the tray cover; the second (upper) door is called the cartridge cover.


Opening the tray cover

- Stops any tray movement immediately.
- Prevents execution of any programmed events which require tray movement.
- Aborts a method when a step requiring tray movement is encountered.

Opening the cartridge cover

- Shuts off the high voltage if it is on.
- Turns off the pump circulating the capillary coolant.
- Moves the detector filter wheel to the closed position.

Figure 5 Cartridge and Tray Covers



1. Tray Cover
2. Cartridge Cover

UV Detector Optics

The UV optics includes an ultraviolet light source, selectable wavelength filters, aperture, capillary, and a single photodiode detector, as shown in Figure 6.


The UV source is a deuterium lamp with a wavelength range of 190 to 600 nm. Two lenses focus and direct the output of the lamp through one of the wavelength-selecting filters located in a rotating wheel behind the capillary cartridge. The beam continues through the aperture in the cartridge plug and through a section of the capillary that has been treated to remove the polyimide coating (the detection window). The non-absorbed beam then continues through a fiber optic cable to a photodiode. The light signal is converted to an electrical signal, digitized, and sent to the P/ACE MDQ computer for processing by the software. This signal is also available as an analog output through a connection on the left side of the instrument.

The design of the instrument insures that the optical system remains in alignment. No user alignments are required.

There are eight positions on the UV filter wheel. UV detector systems are shipped with four standard filters: 200 nm, 214 nm, 254 nm, and 280 nanometers (10 nm bandwidth). These are installed in positions 2, 3, 4 and 5, respectively, of the filter wheel; position 1 is opaque and serves as a shutter for the detection system.

Additional wavelengths are obtained by placing the appropriate filters in positions 6, 7, and 8; the standard filters can also be replaced, if desired. If the instrument will be used with a PDA detector, position 8 should always be left open (no filter). The filter wheel will accommodate widely available ½ inch (12.7 mm) diameter filters with wavelengths from 190 to 600 nm.

Figure 6 UV Optics Layout



901013L.AI

- | | |
|---------------------------|----------------------------------|
| 1. Capillary Aperture | 7. Motor |
| 2. Fused Silica Lenses | 8. 8 Position Filter Wheel |
| 3. Deuterium Lamp | 9. Filter Position (i.e. 214 nm) |
| 4. Lamp Power Supply | 10. Fiber Optic Cable |
| 5. Photodiode | 11. Fiber Optic Connector |
| 6. Fiber Optic Connection | 12. Capillary |

The Photo Diode Array (PDA) Detector

The Photo Diode Array detector, like the UV detector, uses the absorbance of light to detect the presence of samples as they pass through the detection window. Unlike the UV detector, the PDA detector can provide spectral analysis of samples. Spectral signatures obtained in this way can be useful in identifying unknowns.


The PDA detector utilizes the same cartridge configuration as the UV detector. Please refer to the section on the UV detector for a description of the cartridge and to the 32 Karat version 4.0 Installation and Maintenance Manual for instructions on rebuilding the cartridge.

In PDA detection, the full spectrum of light from the deuterium lamp illuminates the capillary (see Figure 12). Light that is not absorbed by samples is delivered by a fiber optic cable to a grating that breaks the light into a spectrum. This spectrum is projected onto an array of 256 photodiodes. With this arrangement, the absorbance profile of the sample can be measured. The PDA detector also allows the simultaneous measurement of light at different discrete wavelengths. By collecting spectral data over time, 3-dimensional plots of absorbance can be created. The photo diode array converts the light signal into an electrical signal. This is digitized and sent to the P/ACE MDQ computer for processing by the software.

The PDA detector always uses filter wheel position #8. It is essential that no filter be in position #8 when the PDA is in use.

The PDA detector is calibrated in wavelength using the discrete emission bands generated by a mercury lamp. The mercury lamp is an integral part of the detector system. Calibration is performed automatically by the software when requested by the user.

Figure 7 Diode Array Optics Layout



- | | |
|--------------------------------|---|
| 1. Capillary Aperture | 11. 9 by 200 μm Fiber Array (Slit) |
| 2. Fused Silica Lenses | 12. Fiber Optic Connector |
| 3. Deuterium | 13. Mercury Calibration Fibers |
| 4. Lamp Power Supply | 14. Mercury Lamp Power Supply |
| 5. Concave Holographic Grating | 15. Mercury Lamp |
| 6. 256 Element Diode Array | 16. Y-fiber Optic Cable |
| 7. Motor | 17. Filter wheel in open position |
| 8. 8 Position Filter Wheel | 18. Fiber Optic Connector |
| 9. Monochromator Entrance Slit | 19. Capillary |
| 10. Fiber Optic Connector | |

Laser Induced Fluorescence Detector

The LIF Detector consists of the LIF Detector Module, the LIF Interconnect Module, and a Laser Module. Refer to Figure 8 for more information. A P/ACE MDQ capillary cartridge with an LIF Detector plug installed is required for use with this system.

The LIF detector uses a laser light source. A 488 nm argon-ion laser and a 635 nm diode laser are available from Beckman Coulter. Other lasers can be adapted. The LIF detector can use dual lasers and dual photodetectors, making it a true dual wavelength system. Refer to the section titled “Using Other Lasers with the LIF Detector” for more information.

A fiber cable transmits excitation light from the laser to the capillary in the cartridge. Substances in the capillary which fluoresce at the laser wavelength are detected. The LIF Detector measures and records this fluorescence, which appears as a peak on the computer screen or printed electropherogram. The LIF Detector can also be used in the indirect mode. In this mode, a fluorescent buffer system is used to detect non-fluorescent components. See Figure 8 and Figure 9.

The initial installation of the LIF Detector will be performed by a Beckman Coulter Field Service Engineer. The P/ACE MDQ can easily be converted between LIF and UV/PDA modes as detector components are modular.

Figure 8 P/ACE MDQ Instrument with LIF Detector






Figure 9 LIF Optical System



The Laser Modules

The following section describes the optional 488 nm and 635 nm laser modules, in terms of how they interface with the P/ACE MDQ instrument and the LIF Detector.


 WARNING	<p>During normal operation of the LIF Detector, laser light is not accessible to the user. To prevent potentially harmful laser light from being emitted from the end of the fiber cable, an interlock mechanism turns off the laser if the laser fiber cable is disconnected from the interconnect module, or if the cartridge cover is opened.</p> <p>Always turn off the laser module and P/ACE MDQ instrument before removing any of the LIF system module(s).</p>
---	--

488 nm Laser Module

The Beckman Coulter 488 nm laser is an air-cooled Argon Ion Laser. Two indicator lights are located on the front of the 488 nm Laser Module. The green light flashes to indicate that the interlock circuit is complete and the laser is preparing to light (in approximately 40 seconds.) Once the laser is operating, the green light will stay on continuously until the laser is turned off, the cartridge cover opened, or until the fiber cable coupler is disconnected from the LIF Interconnect Module.

The yellow light is a current limit indicator. A sensor in the laser module measures the amount of laser light being produced. Circuits in the laser module adjust the amount of current being drawn to maintain a stable 3mW output. The yellow indicator light will come on when the system begins to draw excessive current. This light indicates that the laser tube should be replaced; when it comes on, contact Beckman Coulter Customer Service for assistance.

Figure 10 488 nm Laser Module Back Panel



635 nm Laser Module

The power switch is located on the back panel of the 635 nm laser housing. However, the laser will not light up unless the P/ACE MDQ system is on, and all doors and interlocks are in place. A red light on the front panel of the 635 nm laser module indicates that the laser is on. Refer to Figure 8.

Using Other Lasers with the LIF Detector

Laser units other than the Beckman Coulter Laser Modules 488 and 635 can be used with the LIF detector, as long as the guidelines below are strictly followed.

- If the LIF detector was purchased without a Beckman Coulter Laser Module, an Interlock/Fiber Optic cable (P/N 360671) will have been provided. The custom interlock connector plugs into the LIF Interconnect Module. The fiber optic termination of the cable is the standard SMA 905. The laser unit must be fitted with a fiber launcher which is compatible with this connector.

The two electrical wires in the cable connect to the interlock switches inside the detector and instrument. These switches close when it is safe for the laser to come on. These wires must be connected to the laser's interlock circuit. The connectors and switches in the detector are rated for low voltage (TTL level) signals.



NOTE


To avoid damaging the interlock connector and switches in the LIF Detector Module, do not connect interlock wires in the laser unit to a voltage greater than the rated voltage of the interlock connector and switches (5V, 0.1 A).

U.S. Federal regulations require that this connection be made to limit access to laser energy when the cable is disconnected from the LIF Interconnect Module.

Section 3- Starting 32 Karat Software


To start the 32 Karat Software, select **Programs|32 Karat|32 Karat Software** from the Windows Start Menu button on the Windows Tool Bar.

Figure 11 Start Menu location for 32 Karat Software




32 Karat Software opens to the current location / group window which lists the configured instruments.

Figure 12 32 Karat Software opening window



Tool Bar



The Tool Bar appears at the top and, in some cases, the bottom of the active window. The active buttons within the Tool Bar allow for single mouse click access to many common commands. These buttons are available from the opening window.


Menu Bar

The Menu Bar contains all commands available in 32 Karat Software. The Menu Bar can be accessed by single mouse click or by holding down the Alt key and pressing the underlined letter.

Creating and Configuring an Instrument

In this exercise you will:


- create a connection to a real instrument
- create a virtual instrument that can be used for offline data reprocessing

 NOTE	<p>This section assumes that all hardware and interface boards have been installed. If they have not been installed, please refer to the Installation and Maintenance Manual for instructions.</p>
---	--

Creating a New Instrument


Up to four P/ACE MDQ's can be controlled from a single computer with 32 Karat Software. An unlimited number of virtual instruments can be created for offline method editing or data analysis. Use the following exercise to configure a new instrument. Select the detector type(s) appropriate to your hardware configuration.

1. Right mouse click within the right pane and select **New | Instrument** from the drop down menu.




2. A new instrument will be created with its icon placed in the right pane. Enter the name of the instrument (in this example we will refer to `Instrument 1`) in the highlighted name field. This name will be used throughout 32 Karat Software, in the Instrument window, data reports, and the Instrument logs.

Figure 13 32 Karat Software opening window




Configuring the New Instrument


Configuring the instrument is a process in which the software is prepared to control the hardware and process the data. It consists of 1) identifying the P/ACE MDQ as a unique instrument and 2) identifying the hardware components that are present. Each physical instrument should be represented by at least two icons. One of these will be used for control and data acquisition; the other is used for offline data reprocessing while the first is being used for acquisition. The following example will use a UV detector as a model.

 NOTE	If there is more than one P/ACE MDQ connected to the computer, the configuration process must be repeated for each instrument. Only the instrument that is being configured may be ON during configuration. All other instruments must be turned OFF. Once all instruments are configured, they may all be turned ON.
---	---


1. Right mouse click on the Instrument 1 icon and select **Configure** from the drop down menu.



2. Select **Beckman P/ACE MDQ** from the *Instrument type* drop-down menu; input a name to identify the instrument.




3. Click on the *Configure* button. The Beckman P/ACE MDQ Configuration window is launched with the modules available for configuration listed in the left pane.




Auto Configuration

32 Karat Software can automatically configure all P/ACE MDQ modules in the Instrument. The instrument must be connected to the PC and turned ON.

1. Click the *Auto Configuration* button at the bottom of the window and listen for tray movement within the MDQ. This indicates that the software is communicating with the instrument. The window below will appear during this process.




2. The detector will appear in the configured modules window.



- Right mouse click on the new icon and select *Open*. Verify that the instrument was properly detected and configured.


Figure 14 P/ACE MDQ Instrument configuration for UV Detector



- Each P/ACE MDQ connected to the computer must have a unique device ID. The default setting is 1. If multiple MDQ instruments are installed, select a unique ID (1, 2, 3 or 4) for each instrument. After the device address is selected, click *SET BUS ADDRESS* to send this information to the instrument.
- The type of buffer tray and sample tray present is detected at startup. Tray type is detected each time the tray cover is opened and closed. The instrument autodetects three types of trays: 36 vial buffer trays, large reservoir buffer trays and 48 position sample trays. 96 well plates are not autodetected by the instrument and must be manually configured from this window. Tray configuration may be manually changed at any time.
- When a 96 well plate is selected as the sample tray type, the Enable Tray Definition check box will become active. If this check box is selected, the user has the option to define the height and depth of the plates installed. This feature allows the use of deep-well and other non-standard plates. If this check box is not selected, only standard Beckman Coulter 96 well plates may be used.
- If a Laser Induced Fluorescence Detector has been detected, the LIF Calibration Wizard will become active. This feature is described in greater detail in the Installation and Maintenance Manual.

8. The Filter dialog is used to define the specific filters that are installed in a UV detector. The data must be entered manually. It is important that the values entered match the installed filters exactly. If a PDA detector is installed, position 8 must be empty (no filter installed), indicated by a “0” value.
9. Pressure units may be defined in psi (pounds per square inch) or in mbar (millibars). The units selected will be used only for this instrument.
10. Temperature control indicates the installation of a sample storage unit. If installed, the storage unit may be disabled by selecting “Unavailable” from the drop-down list box.
11. Click *OK* to return to the Instrument Configuration window.

Figure 15 Instrument Configuration screen for the Beckman P/ACE MDQ




Configuration Options

Clicking the Options button will open the Configuration Option dialog. See Figure 16.

Analysis Options

Figure 16 P/ACE MDQ Instrument Configuration Options Dialog box (Analysis Options)



Analysis options define what software functions will be available when data is reprocessed.

PDA allows the analysis of multichannel data from the Photo Diode Array detector.


System Suitability enables automatic review of results. Results outside selected ranges can trigger defined responses.

Qualitative Analysis enables the identification of peaks by migration time, relative migration time or mobility.

Caesar Integration is used to detect peak start and stop. This method is useful for peaks that have abrupt transitions from baseline to peak. It is also useful when S/N (signal noise ratio) is low. Caesar Integration is the preferred method for detecting CE peaks. When un-checked, peak start and stop will be based on slope threshold.

Instrument Options

Figure 17 P/ACE MDQ Instrument Configuration Options Dialog box (Instrument Options)



Standard CE, the time for a peak to reach the window is called Migration Time.

CEC/LC, the time for a peak to reach the window is called Retention Time.


This selection will be reflected on electropherograms and in reports.

12. Click *OK* to exit and save the module configuration.
13. Click *OK* to return to the Opening Menu.



Manual Configuration


Manual Configuration of the UV Detector

Manual configuration may be used to create an instrument when the hardware is offline or in use.


 NOTE	<p>Many of the exercises in this manual require that a virtual UV instrument be configured. Use this procedure to create a UV instrument, regardless of the type of detector you have installed.</p>
---	--

1. From the main 32 Karat Software screen, create a new instrument and open the configuration window as previously described in “Creating a New Instrument” on page 3. Name the instrument “UV for manual”.
2. For this exercise a UV Detector will be configured. Select the UV Detector


 in the left pane and click on the  . The P/ACE MDQ UV Detector icon is added and must now be configured.



- Right mouse click on the new icon and select Open. The P/ACE MDQ Instrument Configuration dialog box is displayed.



- Specify the type of sample trays the instrument is using.
- Click *OK* to exit and save the module configuration.

Configuration Options

For information on configuration options, please refer to “Configuration Options” on page 9.

Manual Configuration of the PDA and LIF Detectors

Manual configuration of the PDA and LIF detectors is performed by selecting the appropriate icons in place of the UV detector. Only one internal detector is allowed per instrument.

Section 4-Direct Control



Introduction


The Direct Control screen is essentially the front panel of the instrument. It allows you to select and change operating parameters with immediate execution. Direct Control is useful in preparing for Operations, Maintenance and Troubleshooting.

The Direct Control graphical interface is used to control the instrument by mouse clicking on “hot” areas of the window. These areas either activate the associated dialog box - allowing modification of the instrument settings - or activate the task directly. Direct control displays an animation of current instrument activities and instrument status in real time.

Direct Control Window

To access the Direct Control window, go to the Control Menu and select **Direct Control | View**.

Figure 18 Hot spots on Direct Control window



1	Capillary Temperature	9	Trays to Load Position
2	Tray Up/Down	10	Injection Dialog
3	Label dialog	11	Stop Current Processes
4	Tray Up/Down	12	Sample Storage Temperature
5	Lamp On/Off	13	Voltage Dialog
6	Detector dialog	14	Current Dialog
7	Autozero	15	Power Dialog
8	Trays to Home Position	16	Pressure Dialog
		17	Vial Position Dialog

The UV Detector is shown in Figure 18 above. The PDA and LIF Detectors are similar. The functions are defined in Table 1 and Table 2.

Table 1 Direct Control Functions


CONTROL FUNCTIONS IN DIRECT CONTROL			
Temperature	Capillary Temperature	Hot Spot	Opens a dialog to set the coolant temperature
	Sample Storage	Hot Spot	Opens a dialog to set the sample storage temperature
Tray Motion	Up/Down	Hot Spot	Toggles trays up and down
	Vial Position	Hot Spot	Opens dialog to move a selected vial to the end of the capillary
	Load	Button	Brings all trays forward for loading
	Home	Button	Returns all trays to their home position
Optics	Detector Control	Hot Spot	Opens dialog appropriate for the configured detector
	Lamp/Laser Status	Hot Spot	Opens dialog allowing lamp or laser(s) to be toggled on/off
	Autozero	Button	Adjusts the detector output to zero
Power Controls	Voltage	Button	Opens dialog to initiate voltage separation
	Current	Button	Opens dialog to initiate current separation
	Power	Button	Opens dialog to initiate power separation
	Pressure	Button	Opens dialog to initiate pressure separation; also used for capillary rinsing
Other Controls	Inject	Button	Opens dialog to initiate injection of sample into the capillary
	Stop	Button	Stops all instrument operations
	Capillary Information	Hot Spot	Opens a dialog to input capillary/ cartridge ID information (used for display only)

Table 2 Status Display

STATUS DISPLAY IN DIRECT CONTROL		
Status	Text	Indicates current instrument function and state
Time Remaining	Text/Graphic	Time remaining for current instrument process
Tray Position	Text/Graphic	Indicates the vial or well at which the capillary end is located
Cartridge Status	Graphic	Indicates if cartridge is installed
Cover Status	Graphic	Indicates if door is opened or closed
Cartridge Temp	Text/Graphic/Hot Spot	Current coolant temperature is displayed
Storage Temp	Text	Current storage temperature is displayed
Storage Status	Graphic	Indicates presence of trays in storage unit
Lamp/Laser Status	Graphic/Hot Spot	Indicates Lamp/Laser(s) on or off
Voltage	Text/Graphic/Button	Indicates present, target, and maximum voltage
Current	Text/Graphic/Button	Indicates present, target, and maximum current
Power	Text/Graphic/Button	Indicates present, target, and maximum power
Pressure	Text/Graphic/Button	Indicates present, target, and maximum pressure
Pressure	Type/Direction Graphic	Indicates positive pressure or vacuum, and capillary end to which the pressure is applied
Polarity	Graphic	+/- indication of electrode potential
Detector signal	Text	Real time signal from installed detector
Capillary Info	Text/Hot Spot	Displays the capillary ID information

Exercise: Conditioning the Capillary

This exercise will use the Direct Control functions to condition a bare fused-silica (BFS) capillary. This capillary will be used for additional exercises in this manual.

 NOTE	<p>This procedure is intended for bare fused-silica capillaries only. Internal capillary coatings may be destroyed by this procedure.</p>
---	---

Materials needed:

- 75 micrometer I.D. capillary, 60 cm total length (50 cm to the detector) installed in cartridge appropriate for detector type configured. See the Installation and Maintenance Manual for instructions on building a new cartridge or installing a capillary.
- Methanol (HPLC grade)
- 0.1N HCl in water
- Capillary Regenerator Solution A (1M Sodium Hydroxide P/N 338424)
- Run Buffer A (P/N 338426)
- Distilled or deionized water (HPLC grade)
- 2 mL vials and red caps
- 36 position buffer trays

Prepare one vial each of Methanol, HCl, Regenerator, and Water. Prepare two vials of Run Buffer A. Prepare one empty vial. Vials should be filled to the shoulder as indicated in the Installation and Maintenance Manual. All vials (including the empty vial) must be capped. Place the vials in the buffer trays as indicated in Table 3.


Table 3 Vial Loading Positions

SOLUTION	Left (inlet) buffer tray	Right (outlet) buffer tray
Methanol	B1	
0.1 N HCl	C1	
Regenerator Solution A	D1	
Water	E1	
Run Buffer A	A1	A1
Empty vial		B1

Check that the capillary cartridge and buffer trays are properly installed (see the Installation and Maintenance Manual). Close the doors. Note that the graphic in Direct Control indicates the presence of the cartridge and the closed door. You should be able to hear the coolant circulating.

Click on the Pressure “hot spot” on the Direct Control screen. The dialog in Figure 19 will open.

Figure 19 Pressure Button dialog



This dialog is used to set up the conditions for a rinse or a pressure separation. It has the following options:

Pressure: Input the desired pressure in psi. The valid range for positive pressure is 0.1 to 100 psi; for vacuum is 0.1 to 5 psi.

Duration: Specifies how long the pressure will be applied.

Direction: Designates whether the flow of fluid will be Forward (from inlet to outlet or left to right) or Reverse (from outlet to inlet or right to left).

Pressure Type: Specifies whether positive pressure or vacuum will be used. This choice defines the valid pressure range.

Tray Positions: Specifies which vials or wells will be located at the capillary ends during the operations. To select the positions, click the Trays button to open the Tray Selection dialog (shown in Figure 20). Click on the appropriate positions to select them. For this exercise, click *B1* on the left tray and *B1* on the right tray to select the Methanol and Empty vials. Click *OK* to return to the Pressure Settings dialog box. The tray positions area should now indicate Inlet:BI:B1 and Outlet:BO:B1.

Type 25 in the Pressure box and 1.0 in the Duration box. Pressure direction should be Forward and Pressure Type should be Pressure. The dialog should now look like the one shown above in Figure 19.

Click *OK*. The trays will move to the designated positions and the rinse will begin automatically. Observe the Direct Control screen during this process.

Perform additional rinse steps according to Table 4. It will only be necessary to change the parameters indicated.

Figure 20 Tray Selection dialog





Table 4 Rinse Step Parameters

Solution	Inlet (L) Position	Outlet (R) Position	Pressure	Time
Methanol	B1	B1	25	1
Water	E1	B1	25	0.5
0.1 N HCl	C1	B1	20	2
Water	E1	B1	20	0.5
Regenerator A	D1	B1	20	2
Water	E1	B1	20	0.5
Buffer A	A1	B1	20	2

At the end of these steps, the capillary has been regenerated and filled with run buffer. The next step is to test the electrical conductivity of the capillary. Click on the Voltage hot spot to open this dialog box:

Figure 21 Voltage Button dialog



This dialog is used to configure a voltage separation.

The parameters are:

Voltage: The voltage at which the separation will be run (0.1 to 30 kV)

Duration: The total time the voltage will be applied. (0.1 to 999.9 minutes)

Ramp Time: How long the power supply will take to reach the voltage (0.10 to 999.9 minutes). Ramp time cannot exceed the Duration time.

Voltage Max: The maximum voltage the system will be allowed to reach. This must be greater than or equal to the set voltage.

Current Max: Maximum current the system is allowed to reach. If this value is reached at a voltage lower than the set voltage, the set voltage will not be achieved. Maximum allowed current is 300 μ A.

Tray Positions: This is the same as previously described under Pressure.

With Pressure: Allows the simultaneous application of pressure and voltage. When selected, additional parameters of pressure level and direction become available. Maximum pressure is 100 psi. Pressure can be applied to the inlet (Forward), outlet (Reverse), or both ends of the capillary.

With vacuum: Allows the simultaneous application of vacuum and voltage. When selected, additional parameters of vacuum level and direction become available. Maximum vacuum is 5 psi. Vacuum can be applied to the inlet (Forward), or outlet (Reverse), end of the capillary, but not both.

External Adapter: This box must be selected when the external adapter accessory is in use. For additional information, see the Installation and Maintenance Manual.

Polarity: Determines the sign of the charge on the electrodes. The graphic will indicate the selection.

For this exercise, use the Trays button and select A1 on the inlet side and A1 on the outlet side. Both ends of the capillary will be in Run Buffer A. Set the other parameters as follows:

Voltage: 30kV

Duration: 2 minutes

Ramp time: 0.5 minutes

Click *OK*. The trays will move to position the selected vials at the ends of the capillary. The voltage will begin to ramp up to 30 kV. Observe the Direct Control screen during this process, taking particular note of the Voltage, Current, and Power displays. The Current should reach between 27 and 33 mA when the Voltage is at 30 kV. At the end of two minutes the voltage will return to zero and the Status will return to idle.

This completes the Direct Control exercise. It may be helpful to experiment with the other Direct Control functions before proceeding to the next section.

Section 5-Method Editing

Creating a Method

Introduction

In the section on Direct Control the exercise on regenerating a capillary required user input at each step of the process. A *method* combines a series of steps into a logical process. The steps in a method are executed automatically by the instrument. This section will describe the basic steps in creating an instrument method for data analysis. There are many more features, such as data analysis and reporting, that can be incorporated into a method. Some of these are covered later in this manual. For more information, see Appendix 1.

To create or edit a method you must be in the Instrument Window. To open this window, go to the 32 Karat main screen. Right click on the icon you created to represent your instrument system and select *Open Offline*. The instrument window will open after a few seconds. When the window opens, the Instrument Wizard will appear. It is possible to open the method editing dialogs by selecting *Create or Modify a Method* from this dialog. For now, just click *OK* to close the Instrument Wizard.

We will write a new method to run an instrument test sample. To begin, select **File | Method | New** from the menu bar. The name of the method in the Instrument Window title bar will change to “untitled.met.” To access the instrument control and data acquisition sections of the method, select **Method | Instrument Setup** from the menu bar. A window will open that contains three or four tabs (depending on your detector type). Select the tab marked “Initial Conditions” to bring it to the front. The dialog in Figure 22 will appear. This dialog is used to set instrument parameters at the start of a method, before the separation process begins.

Initial Conditions Tab


In the previous section, **Direct Control**, the exercise on regenerating a capillary required user input at each step of the process. This section shows how to combine a series of these steps into a logical process called a *Method*. The steps in a method are executed automatically by the instrument. This section will describe the basic steps in creating an instrument method for data acquisition. There are many more features, such as data analysis and reporting, that can be incorporated into a method. Some of these are covered later in this manual. For more information, see Appendix 3.

In this exercise you will write a method to run an instrument test mixture. To write a new method, select **File | Method | New** from the menu bar.

To access the instrument control and data acquisition parts of the method, select **Method | Instrument Setup** from the menu bar. A window will open that contains

three or four tabs, depending on the detector configuration. Select the tab marked *Initial Conditions*. The dialog in Figure 22 will appear.

Figure 22 Instrument Setup dialog - Initial Conditions tab



Auxiliary data channels

The 32 Karat Software gives you the option to collect any one or all instrument parameters: voltage, current, power, and/or pressure. Select the channels you wish to record by clicking the appropriate box(es).

The max kV and max μA boxes are used to set the allowable limits for these parameters. Voltage and current are interrelated by the expression $V=IR$. The system will limit both parameters whenever one limit is reached. For example: assume a voltage set at 30kV and a current limit set at 10 μA . With some buffer systems, a voltage of 12 kV will generate a current of 10 μA . In this case the voltage will not exceed 12 kV, as the current limit will be the determining factor.

Mobility channels

This topic is covered later in this manual.

Temperature

Sets the initial temperature of the cartridge coolant and the sample storage unit, if installed. These settings can be changed by the time program.

Peak detect parameters

These parameters are used to trigger specific events such as fraction collection. They are not used in peak integration or data analysis; peak integration is discussed later in this manual.

Trigger settings


The P/ACE MDQ can be forced to wait until certain conditions are met before beginning a run. These are selected by clicking the appropriate box. If “Wait for external trigger” is selected, the MDQ becomes a slave device, and will not start until an external signal is received. The two “Wait for temperature” options assure that the system has reached the correct operating temperature before beginning a run. These options delay only the start of the time program. Parameters set in initial conditions will occur without a wait.

Inlet and Outlet trays

The types of trays that will be used when a method is run must be specified here. When the method is run, this information will be compared to the tray types that are configured in the instrument. If there is a tray type mismatch, the method will not run.

Detector Initial Conditions


For this exercise, set conditions in the Initial Conditions tab to match Figure 22.



UV Detector Initial Conditions

Click on the UV Detector Initial Conditions tab.

Figure 23 Instrument Setup dialog - UV Detector Initial Conditions tab



Electropherogram channel

The Acquisition enabled box must be selected in order for data to be collected and saved.

The Wavelength box allows the selection of one of the filters installed in the filter wheel. The available selections will be those set up during instrument configuration.

Data rate is selectable from 0.5 to 32 Hz (data points per second). Narrow peaks require a higher data rate. Data rate interacts with the Peak width parameter, set under Filter settings in this dialog box.

Filter

This filter refers to an algorithm for filtering the collected data, and does not refer to the wavelength filters. Data filtering is necessary to remove extraneous noise that might interfere with data analysis. Filter setting requires two steps. The selection of the filter type (High Sensitivity, Normal, or High Resolution) and the setting of the Peak width. High Sensitivity increases signal to noise at the expense of resolution. High Resolution increases resolution at the expense of signal to noise. Normal is a compromise setting that is suitable for most analyses. The filtering algorithm is most efficient if it has information about the number of data points that make up a peak. The optimum number of points is from 16-25. There are also options for “less than 16” and “more than 25” data points per peak. The Data rate setting in the Electropherogram channel dialog can be adjusted so that a typical peak has from 16-25 data points. The other Peak width settings should be used only if changing the Data rate is not desirable.



NOTE

Peak width here is not the same as Peak Width in data integration.

Relay 1 and Relay 2

The two relays incorporated into the P/ACE MDQ can be used to control or signal other devices. This dialog allows you to set the state of the relays at the start of the method to “Off” (open) or “On” (closed). These relays may also be controlled from the time program. See the Installation and Maintenance Manual for more information on using the relays.


Absorbance Signal

This parameter has two options, Direct and Indirect. Direct is used when the analytes have higher UV absorbance than does the background electrolyte. In this case, upward deflecting peaks are produced when analytes pass through the detector window. In some modes of CE, the background electrolyte has a higher absorbance than the analytes. In this type of separation, the passage of analytes creates negative peaks (valleys) in the electropherogram. By selecting Indirect, the entire electropherogram is inverted, so that these negative peaks appear as positive peaks.

For this exercise, set up the UV Detector Initial Conditions dialog to match Figure 23.

PDA Detector Initial Conditions

Figure 24 Instrument Setup dialog - PDA Detector Initial Conditions tab



Electropherogram Scan Data

This dialog controls the parameters for the collection of 3D data. The Acquisition enabled box must be selected in order for data to be collected and saved. Data rate is selectable from 0.5 to 32 Hz (data points per second). Narrow peaks require a higher data rate. Data rate interacts with the Peak width parameter in Filter settings (below). Scan data collected at high data rates can result in very large data files. Scan range specifies the wavelength range over which data will be acquired (limit 190-600 nm).

Filter

See the UV Detector Initial Conditions section for information on this topic.

Relay 1 and Relay 2

See the UV Detector Initial Conditions section for information on this topic.

Reference Channel

See the sources in Appendix 1 for information on this feature.

Absorbance Signal

See the UV Detector Initial Conditions section for information on this topic.

Electropherogram Channel Data

Up to three independent data channels can be collected simultaneously with the scan data. For each channel the following must be input:

Acquisition enabled: Check this box to enable the collection of data in this channel.

Reference channel: See the sources in Appendix 1 for information on this feature.

Wavelength: This parameter specifies the central wavelength for the data channel.


Bandwidth: This parameter specifies the width, in nanometers, of the data collection channel. Larger values will give higher signal to noise ratios than will low values. Low values will be more specific for molecular structure than will higher values.

Peak detect: See the sources in Appendix 1 for information on this feature.

For this exercise, setup the PDA Detector Initial Conditions to match Figure 24.

Multi-Egram


Figure 25 Instrument Setup dialog - Multi-Egram tab



This tab allows for the collection of additional data channels which are extracted from the 3D data file. See the sources in Appendix 1 for information on this feature. It will not be used in this exercise. The Multi-Egram tab is only present when a PDA detector is installed.

LIF Detector Initial Conditions

Figure 26 Instrument Setup dialog - LIF Detector Initial Conditions tab



Electropherogram Channel 1 and 2

The LIF Detector is available in one and two channel versions. The programming of these channels is identical. In this exercise we will use only one channel. If your system is configured with non-Beckman Coulter, Inc. lasers, the test mixture used in this exercise may not perform as described.

Acquisition enabled is checked to enable the collection of data in this channel.

Dynamic range specifies the upper limit of the signal range. Higher values allow for the collection of larger peaks without truncation, but at the cost of sensitivity. Lower values give greater sensitivity, but large peaks may be truncated. This value should be slightly larger than the largest peak expected.

Signal

Select Direct when separating fluorescence samples (peaks) in a non-fluorescent electrolyte background electrolyte. If a fluorescent background electrolyte is used, non-fluorescent samples may be detected. Selecting Indirect in this case will invert the signal so that negative peaks appear as positive peaks.

Data Rate

Specifies the rate at which data will be acquired. Both LIF data channels will have the same data rate. They cannot be specified separately.


Relay 1 and 2

See the UV Detector Initial Conditions section for information on this topic.

For this exercise, setup the LIF Detector Initial Conditions to match Figure 26.

Time Program

Figure 27 Instrument Setup dialog - Time Program tab



The Time Program window is arranged like a spreadsheet. An event is entered into each line. Events will be executed in order, top to bottom. Each line in the window contains these columns:

Time: the point after time zero at which the event will occur.

Event: the action that will occur (see below).

Value: this will vary depending on the action selected.

Duration: time the event will last.

Inlet and Outlet Vial: where the capillary ends will be during the event.

Summary: a system generated description of the event.

Comments: a user generated annotation of the event.

Time is not a required event. Events that have no time associated with them will be run in the listed order, top to bottom, and each will be finished before the next event begins. Timed events must be grouped together; a group of timed events cannot be interrupted by an untimed event. Untimed events can only occur before a group of timed events. Some events do not have a time option, others can be timed or untimed. Data acquisition begins with the first timed event (time 0.00); it ends when the method ends or when a STOP DATA event is reached.



To program a line, click in the EVENT box. Clicking on  opens a menu of events. Selecting an event will open a dialog box for that event.

Figure 28 Time Program Event List




The following will describe only those events that are used in this exercise. For information on the other event types, see the resources in Appendix 1.

Separate

The separate dialog is used to control the conditions under which the separation process takes place. Every method must have at least one separation step. The step at time = 0.00 is usually a separation.

Figure 29 Separate dialog



It has the following parameters:

Separation Type: Electrically driven separations can be done at controlled Voltage, Current, or Power. When one of these is selected, the other two will float to a value determined by the resistance of the capillary contents. Voltage and Current cannot exceed the limits sets in the Initial Conditions window. Separations can also be programmed to use pressure or vacuum to move the fluid in bulk through the capillary. Voltage, Current, or Power can be combined with Pressure or Vacuum so that two processes are at work simultaneously.


Polarity determines the direction of the current. The charge on the electrodes is indicated by the graphic in this dialog.

Values allows for input of the set points for the Separation Type parameters. The available options will change depending on the Separation Type selected. Ramp time is only valid for electrical separations. It determines how long it will take for the Voltage, Current or Power to change from the present to the programmed level.

Tray Positions can be selected graphically by clicking on the Trays button. The type of trays shown are determined by the settings in the Initial Conditions dialog. When a method is used in a Sequence Table it may be desirable to change the vial positions after a specified number of runs. The inlet, outlet, or both positions can be incremented

automatically by selecting the appropriate boxes and inputting the desired number of runs between changes.

Figure 30 Tray Selection dialog



Pressure Direction specifies whether pressure or vacuum will be applied to the inlet, or outlet end of the capillary. Pressure (but not vacuum) can also be applied to both ends of the capillary at the same time.


At Time specifies that this will be a timed event at the time specified. Separation is usually a timed event.

External Adapter changes the way the instrument manages the power supply. Select this box *only* if the External Adapter accessory is in use. See the Installation and Maintenance Manual for more information.

Rinse

The rinse event is used to clean the capillary and to load fresh buffer or other separation media.

Figure 31 Rinse dialog



It has the following parameters:

Pressure Type selects the mechanism to be used to move fluid through the capillary.

Tray Positions functions exactly as described in the Separation event.

Values specifies the magnitude of pressure to be delivered, and for how long.


Pressure Direction specifies whether pressure will be applied to the inlet or the outlet end of the capillary.

At Time functions exactly as described in the Separation event.

Inject

The inject event is used to deliver a precisely measured amount of sample into the capillary. This step is always untimed, and usually precedes the first separation step.

Figure 32 Inject dialog



It has the following parameters:

Injection Type: Sample can be delivered to the capillary by positive pressure, vacuum or by the application of voltage (electrokinetic injection).

Polarity specifies the charge on the electrodes during a voltage injection.

Pressure Direction specifies whether pressure or vacuum will be applied to the inlet or the outlet end of the capillary.

Values specifies the magnitude of the pressure, vacuum, or voltage, and how long it will be applied. Higher values will inject more sample.


Tray Positions functions exactly as described in the Separation event.

Sequence Table: When a method is used in a Sequence Table, certain parameters in the inject event can be overridden by values input in the Sequence. This parameter determines whether the method Inject event or the Sequence Table has priority.

Multiple injection events are allowed. If multiple injections are used, only the first injection event may be overridden in the sequence table or incremented

Autozero


Figure 33 Autozero dialog



This event causes the detector output to be reset to zero. It may be timed or untimed.

Stop Data


Figure 34 Stop Data dialog



Data collection starts at time = 0.00. Data collection will continue until the end of the method unless a STOP DATA event is encountered. This event can be used to avoid data collection during steps such as post-run capillary cleaning. STOP DATA is always a timed event.

End

Figure 35 End dialog



End is an optional event. The method will not continue beyond an end event. It is always a timed event.


In the next part of this exercise we will create a time program to go with the initial conditions previously entered. The example method will do the following:

- Perform a mini-regeneration on the capillary
- Fill the capillary with buffer
- Inject test mix
- Separate the test mix
- Perform a post-run rinse

This method will use a vial tray setup similar to that used in the section on Direct Control. The vial positions are specified in the Table 5.

Table 5

Solution	Left (inlet) buffer tray	Right (outlet) buffer tray
Test Mix	C1	
Regenerator Solution A	D1	
Water	E1	
Run Buffer A	A1, A2	A1
Empty vial		B1

Click on the Time Program tab to begin. The tab opens with a blank spreadsheet. Click in the Event box and select the arrow, . From the drop down menu select Rinse. Click on the Trays button. Select D1 on the Inlet (left) side and B1 on the Outlet (right) side. Before clicking *OK*, right click on one of the buffer positions you just selected. A dialog box will open that allows you to identify the contents of that position. Insert the appropriate comments from the table above. Repeat for the other vials. When done, click *OK* to return to the Rinse window. For this step, accept the default values for all other parameters. Click *OK* to return to the Time Program spreadsheet. The key parameters have been entered into the table automatically. This rinse is an untimed step, so there is no value in the Time window.

In the blank line at the bottom of the table, click in Event column again and select another rinse event. This time, program a rinse with water from position BI:E1 to position BO:B1. Change the rinse time to 0.5 minutes. This event will be added to the table when *OK* is clicked in the dialog box.

Click in Event column again, and add a step to fill the capillary with buffer rinse from position BI:A1 to BO:B1. This time, use 25 psi for 1 minute. Click *OK* when done.

Click in the Event column in the bottom row again. Select Inject as the event type. Select the Trays button, and program the injection from position BI:C1 to BO:A1. Click *OK*. Input 0.4 psi as the pressure and 3 seconds as the time. Click *OK*.

At this point in the method we are ready to separate the components in the test mix. Select a new Event box, and choose Separate as the event type. Use the Trays button to select positions BI:A2 and BO:A1. Click *OK*. Select a Voltage separation. Set the

voltage to 30 kV, ramp time to 0.2 minutes, and duration to 6 minutes. This will be a timed step, so the At time button must be checked. Set the time to 0.00 minutes. Click *OK* to put this event into the table.

We will be using a final rinse in this method, so we need to insert a Stop Data event. Select the new Event box and choose Stop Data as the event type. Enter a value of 6.00 minutes to force data collection to stop at the end of the separation step.

Our last step will be a rinse with water from the water vial (EI) to the empty vial (BI). Program this rinse for 1 minute at 20 psi. Select At Time and use a value of 6.00 minutes.

It is a good practice to Autozero the detector some time after the run starts. In the new event box select Autozero. Select At Time and input a value of 1.00. When this item is entered into the Time Program spreadsheet, it will automatically move to the correct time sequence.

At this point your spreadsheet should look like that in Figure 36.

Figure 36 Instrument Setup - Time Program

	Time (min)	Event	Value	Duration	Inlet vial	Outlet vial	Summary	Comments
1		Rinse - Pressure	20.0 psi	1.00 min	BI:D1	BO:B1	forward	
2		Rinse - Pressure	20.0 psi	0.50 min	BI:E1	BO:B1	forward	
3		Rinse - Pressure	25.0 psi	1.00 min	BI:A1	BO:B1	forward	
4		Inject - Pressure	0.5 psi	4.0 sec	BI:C1	BO:A1	Override, forward	
5	0.00	Separate - Voltage	30.0 KV	6.00 min	BI:A2	BO:A1	0.20 Min ramp, normal polarity	
6	1.00	Autozero						
7	6.00	Stop data						
8	6.00	Rinse - Pressure	20.0 psi	1.00 min	BI:E1	BO:B1	forward	
9								

The method can still be edited at this point. Items in the Time, Value, Duration, Inlet, and Outlet columns can be edited by clicking on the existing value and typing in a new value. Alternatively, clicking on the name of any existing event will open the drop down menu of event types. Select Current Event to re-open the event dialog, or select a new event to change the event type.

It is very important to **SAVE YOUR METHOD**. Select **FILE | METHOD | SAVE AS** from the menu bar. Give the method a name which will be meaningful. For example:

My Test Method of 5-16-2000

Click Save. The default path at installation is:

d:\32Karat\Projects\Default\Methods.

Your system administrator may have assigned you to a different default path.

The method is now ready to run.

In the next section you will run the method you just created. For now, practice writing and editing methods to explore the various parameters and how they interact. The resources in Appendix 1 describe the events and parameters not covered here.

Section 6-Running the System

Single Run

In the last section you created an instrument method. In this exercise you will run that method in Single Run mode. In this mode, the method must be manually started by the user prior to every run. Single mode is useful in method development, where the results of one run will suggest modifications to the method or other procedures. In the second part of this section you will create a Sequence to do multiple runs.

You must be “online” to run samples and acquire data. If you are not online, “(offline)” will appear in the title bar of the instrument window. Close any offline windows before proceeding. From the main 32 Karat Software screen, double click on an instrument icon to go “online”.



NOTE

The detector lamp or lasers must be ON prior to starting a method. The status of the lamp or lasers is available from Direct Control. The UV lamp may be started there. Lasers may be started from Direct Control if the power switch on the laser module is already ON and all safety interlocks are closed.

Materials needed:

- Cartridge prepared in Direct Control exercise (75 micrometer id capillary, 60 cm total length (50 cm to the detector) installed in cartridge appropriate for detector type configured. (See the Installation and Maintenance Manual for instructions on building a new cartridge or installing a capillary.)
- Capillary Regenerator Solution A (1M Sodium Hydroxide P/N 338424)
- Run Buffer A (P/N 338426)
- Distilled or deionized water (HPLC grade)
- 2 mL vials and red caps
- 36 position buffer trays
- For UV or PDA Detection: Test Mix B (P/N 501333)
- For LIF detection using a 488 nm argon-ion laser:
 - LIF Detector Test Mix (P/N 477615)
- For LIF detection using a 635 nm diode laser:
 - LIF Detector Test Mix (P/N 726022)

Load the trays as follows:

Table 6

	Left (inlet) buffer tray	Right (outlet) buffer tray
Test Mix*	B1	
Regenerator Solution A	D1	
Water	E1	
Run Buffer A	A1, A2	A1
Empty vial		B1
* P/N 501333 for the UV or PDA Detector P/N 477615 for an LIF Detector with a 488 nm laser P/N 726022 for an LIF Detector with a 635 nm laser		

Install the loaded trays and the capillary cartridge into the instrument as described in the Installation and Maintenance Manual. If the instrument software window is not already open, open it as described in the section **Starting the Software**.



Select **Control | Single Run** from the menu bar or click on the Single Run icon  to open the Single Run dialog shown below.

Figure 37 Single Run Acquisition dialog



You must provide at least a: Sample ID, Method, and Data File. The Data File will be stored in the directory identified in the Data Path. Sample ID can be any text string that identifies your sample. This ID will appear on all subsequent reports. The method file can be selected by clicking on the Open icon and navigating to the desired method. For this exercise, select the method you created in the Method Editing section.

Data path is selected in the same way as the selection of a method. In this case, you are selecting a directory rather than a file. The data file is most easily entered by typing a unique name in the box. The file name must be a name that does not already exist in the data directory; if the name already exists, the run will not start.

The remaining items on this screen are not required to run a method. For information on those features, see the resources in Appendix 1.

When ready to begin, press Start. The system will examine the method to verify that it is appropriate for the instrument configuration. The method is then downloaded to the P/ACE MDQ. The instrument will first perform a brief start-up check, and then will execute the method. This is a good opportunity to become familiar with how the instrument operates. Knowledge of normal operations is essential when things go wrong. There is a possibility that your first run may not perform as expected. If error messages appear, or the data does not look like the example shown here, check the following.

- Are the vials correctly filled with the proper fluid?
- Are the vials in the positions specified by the method?
- Was the correct method selected?
- Was the method written for the current instrument configuration (e.g., a PDA method cannot run on a UV or LIF system)?
- Is the method correct (compare to Section 5)?
- Do the tray types in the method Initial Conditions match the tray types configured for the instrument?


This method will be used for additional exercises. Please do not continue to the next section until you have achieved a successful single run.

Programming a Sequence

A Sequence is a list of methods and data files that will be used to run a batch of samples without user intervention. Sequences can be used to acquire data (run the instrument) or to batch reprocess existing data files. In this exercise, you will create a Sequence to acquire data from multiple runs of the test mix method. Later in this manual you will use a sequence to perform batch reprocessing of data. You must be “Online” for this exercise.

To create a new sequence, select **File | Sequence | New** from the Instrument window. This will open the Sequence Wizard dialog shown here. The Wizard consists of 5 screens. Not all of them will be used for every sequence, and not every feature of each screen will be used in this exercise. For information on features not covered here, see the resources in Appendix 1.

Figure 38 Sequence Wizard - Method dialog






The first screen requires that you select a method. Click on the  icon and navigate to the method file created in the last section. Under Data File Type select For Acquisition. Amount Values will not be used in this exercise. Click Next when ready.

Figure 39 Sequence Wizard - Unknowns dialog




This screen is used to create the Sample ID and Data File ID. Type a text string in the Sample ID box. The blue arrow to the right of the text box opens the following menu:

Figure 40 Sample ID text options menu




Selecting an item from this menu will cause a symbol to be inserted into the Sample ID field. In the example shown, Line Number was selected. This parameter will be

automatically incorporated into the Sample ID when the sequence is run. You may select any combination of these items.

Data Path specifies the directory where the data files will be stored. Select this by clicking on the  icon and navigating to the desired folder.

Type a text string into the Data File box. The blue arrow to the right of the text box opens this menu:

Figure 41 Data File text options menu




As described above for the Sample ID, this menu is used to insert codes that will be automatically inserted into the filename when the sequence is run. In this example, Date and Time has been selected. (The Open File choice is used to reprocess existing data files, not for data acquisition).

Number of Unknown Runs in Sequence will determine how many lines there are in the Sequence Table. For this exercise, input 3.

Repetitions per Run will determine how many times each line in the Sequence Table will be run. An identifier will automatically be added to the file name for each repetition when this option is used. For this example, input 1. Click Next when ready.

Figure 42 Sequence Wizard - Vials dialog



If the method being used allows vial position override (see the Method exercise), the desired starting position can be input here. If override is not allowed or if auto-incrementing has been set, any input here is ignored. A set of calibration vials can also be identified in this window. See the resources in Appendix 1 for more information on calibrating through the Sequence Table.

Figure 43 Sequence Wizard - Calibration dialog

Sequence Wizard - Calibration

Calibration ID :

Calibration path :

Calibration file :

Number of calibration levels :

Repetitions per level :

Clear all calibration at start of sequence

Create a separate row in the sequence for each repetition

Multiple calibration sets

Number of unknown runs between sets :


Intersperse calibration vials with unknown vials

Reuse calibration vials from first calibration set

? Cancel < Back Next > Finish

In this example we will not be doing calibration. Calibration is covered in the resources listed in Appendix 1. Input 0 in the Number of calibration levels box, and click Next to continue.

Figure 44 Sequence Wizard - Reports dialog



This screen allows you to have reports generated automatically at the end of the sequence run. We will not generate a report in this exercise. All the other choices should be blank. Click Finish. A window similar to this one will appear.

Figure 45 Sequence Table

Run #	Status	Run Type	Level	Conc Override	Reps	Sample Inject Inlet	Sample Inject Outlet	\$
1		Unknown	0	n/a	1	BI:A1	BO:A1	1.
2		Unknown	0	n/a	1	BI:A1	BO:A1	1.
3		Unknown	0	n/a	1	BI:A1	BO:A1	1.
4								

This screen contains many more columns than can be displayed at one time. Most of this information can be edited. Scroll left and right to review the entire screen. Not all of the columns will be used in this exercise. The important ones are described below.

Run Type - Because we selected to have a Summary Report generated, all of the runs are of the Summary type. Begin and End marks the limits of the runs to be included in the report.

Reps tells the system how many times to run the line. Click in the Reps box for the second line and change the Reps to 2. Line two will now execute twice, giving a total of four runs in this three line sequence.

Scroll to the right until the window resembles the one shown here.

Figure 46 Sequence Table (cont'd)

Run #	Sample ID	Method	Filename	Sample Amt.	ISTD Amt.	Multiplier	
1	Test Mix 001	My Test Method of 5-16-2000.met	Test Runs <D>.DAT	1	1	1	
2	Test Mix 002	My Test Method of 5-16-2000.met	Test Runs <D>.DAT	1	1	1	
3	Test Mix 003	My Test Method of 5-16-2000.met	Test Runs <D>.DAT	1	1	1	
4							
5							

Sample ID is based on the Sample ID text you entered in the Wizard. The line number has been added to the Sample ID you input. Click in the Sample ID box for Run #2 and change it to read Test Mix 002 Twice.

Method shows the method you entered in the Sequence Wizard. To change the method, click in the Method box and click on the green icon. (For this exercise, do not change the method.)

Filename is based on the File Name Text you entered in the Sequence Wizard. The filename still shows the Date and Time symbol <D>, because date and time are unknown until the sequence is run.

When you are finished editing the table, select **File | Sequence | Save As** from the menu bar. Use a Sequence name that has significance to you, for example:

Test Sequence 5-16-2000.

Running the Sequence

You will now run the sequence. The instrument must be prepared for operation as described under Single Run Mode earlier in this section. Verify through Direct Control that the instrument UV lamp or laser is turned on. The test sample and buffer vials must be loaded and installed. You must be “online” before proceeding.




Click the double green arrow  icon from the Instrument window toolbar. This icon will not be available if you are in Offline mode.) The dialog shown here will be displayed:

Figure 47 Run Sequence dialog



Click on the Open  icon and select the sequence you just created. The sequence name will appear in the text box.

Run Range allows you to run only part of the Sequence. For this exercise select All.

Printing allows the system to automatically print reports when they are ready. If you do not have a printer connected, leave this box unchecked.

When complete, the dialog should look like the one shown in Figure 47, but with your sequence name in place of the one shown. Verify that the instrument is loaded and ready to run, then Click Start.

The software will validate that the methods in the sequence are appropriate for the current instrument configuration. If any problems are detected, a message will be displayed and the sequence will not be run. Correct any problems noted and restart the sequence.

After a successful sequence validation, the method in the first line will be downloaded to the P/ACE MDQ and the run will begin. Observe the operation of the instrument as

a guide to future operations and troubleshooting. You may open the Direct Control window during data acquisition to view real-time information on instrument status. During the run, the data will be displayed in real-time in the instrument windows. One window will show the Absorbance (or Fluorescence) signal, and a second window will show the Current signal, as those channels were selected in the method. If you are using a PDA detector, a third window will display a contour plot of the absorbance signal.

At the end of each run, the method for the next run will be downloaded afresh before the new run begins. This feature allows you to make changes to a method while a sequence is processing. The version of the method that is current at download (last saved) is the version that will run. You may also edit the Sequence Table during processing. Existing lines that have not been started may be edited or deleted, and additional lines may be added to the table. (Lines that have already completed or are currently executing can not be edited).

When the sequence is complete, you will have acquired data from four runs. The two replicate runs from line two will have file names to which “rep 1” and “rep 2” have been added.


It will be helpful to experiment with the Sequence Table before proceeding to the next section. Make changes, process the sequence, and observe the effects. You may wish to combine these experiments with the writing of new methods. A series of different methods can be entered into the Sequence Table. When processed, the results will demonstrate the effects of changes in parameters. For example, you might create a series of methods at different separation temperatures (e.g., 20, 22, 24, 26, 28, and 30°C) to examine the effect of that parameter on your separation. Another series of experiments might examine the effects of changing the volume injected. If your method allows override, this can be done with a single method by editing the Sample Inject Duration column in the Sequence Table.

Section 7-Integration

Description

An electropherogram is a graphical depiction of the passage of molecules through the detector window. The area under the peak is directly related to the amount of the analyte present in the sample. The process of determining this area is known as Integration. In practice, integration includes steps that determine the points at which peaks start and stop, and the shape and position of the baseline. If sample components are not fully separated from one another (resolved), it is also necessary to estimate where to divide the area of the two peaks. A very simple data set is shown in the figure below.

Figure 48 Typical data set




901023L.A1

In this example, there are four components. Two are well resolved from the others. Two other components are only partially resolved. The four components differ in peak area. This may indicate that they are present at different levels, that they produce different responses in the detector, or both. Quantitation will be covered in the next section. This exercise will focus on peak detection and baseline construction.

Integrating Data

There are two required parameters for peak integration: width and threshold. More complex separations will require the use of more parameters. The test mix separations used in the previous exercises yield very simple electropherograms. In this exercise we will use examples that are installed with the software. The principles of integration are the same regardless of detector type. Configure a UV instrument for offline analysis.


From the Instrument window select **File | Method | New**. Close the instrument setup window when it appears, as it will not be needed for this exercise.

 NOTE	Review the instructions in the Preface on page ii before using the sample files. Work from a copy of the data files and not the originals.
---	--

From the Instrument window, select **File | Data | Open**. Navigate to your *copy* of the Data Samples folder and select the file named Data Sample 1.dat. The data will open in the Instrument window. Click the electropherogram window to select it, then type `Ctrl-Z` to zoom the data to maximum size. Only Channel A data will be used.


Auxiliary data (current, voltage, power, pressure) can not be integrated. Right click on the electropherogram and select **Annotations**. The following dialog will appear:

Figure 49 Trace Annotations dialog



This dialog allows the selection of items that will be displayed in the electropherogram window. Select the desired items (Peak#, Area, and Migration Time) by clicking the item in the left box, then using the green arrow to move the item to the right box. Verify that Baseline and Show undetected peaks are selected (checked). When the dialog looks like the one above, click OK. (Some items may not appear on the electropherogram display until the data is integrated.)

From the Method menu select Integration Events. This will open the Integration Events spreadsheet. There are two items in this table by default: Width = 0.2 and Threshold = 50. Every Integration Table must contain at least one Width and one Threshold value.

To initially integrate the data, use only the default parameters. Select **Analysis | Analyze** from the menu bar, or click on the Analyze icon  on the button bar. The integration process takes a few seconds; larger data files, and files with more peaks, will take more time to process. During integration the Analyze button will turn into a “Stop” button. Clicking the Stop button will abort the analysis.

When integration is complete, approximately 26 peaks will appear. Above or below each peak the peak number, area, and migration time will be displayed. The default parameters will have identified major peaks, but also captured a number of baseline artifacts that are not of interest. (To zoom in on the baseline to view integration results in detail, you can use the mouse to “rubber box” the display. Click and hold anywhere in the image. Drag the mouse to highlight the desired area. When the mouse button is released, the selected area will fill the window. Use `Ctrl-Z` to zoom out.) Examining the area between 3.0 and 3.4 minutes will demonstrate that even negative peaks (valleys) are being detected.

Optimizing Integration

By selecting **Window | Tile Horizontally** the electropherogram and the Integration events table can be displayed together. Resize the windows as desired by clicking and dragging a window's edge with the mouse. Click in the electropherogram window, then use `Ctrl-Z` to zoom out, if necessary. Click in the Value column of the Threshold event in the Integration Events Table. Enter the number 100. Analyze the data again. This time there should be about 10 peaks visible. Close examination will reveal that these are most likely real peaks, and not baseline artifacts, although some of them are quite small. Examine the area between 14 and 15 minutes, for example.

Threshold determines how high a peak must rise above the background noise before it is recognized as a peak by the integration software. Increasing the value of Threshold removes background artifacts. If taken to extremes, it can also remove real data. Set the Threshold value to 1000 and analyze the data again. You should observe that the minor peaks are no longer integrated, and the baseline below the large peaks has shifted dramatically.

Before examining the effects of the Width setting, restore the Threshold value to 100. Click Analyze again. Zoom in on the base of the major peaks. Change the Width value to 0.1 and Analyze. The tail of the second major peak should now be split off as a

separate area. Change the Width to 0.5 and analyze. The tail is again cut off, but in a different place than it was with the smaller value for the window.

The question will inevitably arise, which is the correct integration? The answer is, it depends. You, the analyst, must use your knowledge and judgement to determine the correctness of the integration results. In this example, the peak tailing is a separation artifact, and the area of the tails should be included in the area of the peak. To obtain this result, with THIS data file, a Width of 0.2 and a Threshold of 100 seem to be appropriate.

Integration Parameters

There are many other integration parameters besides the required Width and Threshold that can be used to generate an integration result that properly reflects the components of the sample. These are discussed at length in the resources listed in Appendix 1. In the next part of this exercise, we will use some of these additional parameters to aid in the integration of a somewhat more complex data file.

To begin, return the Width value to 0.2 and the Threshold value to 50 (the default values). Open the data file Data Sample 2.dat. Use `Ctrl-Z`, if necessary, to view the full file. This data file differs from the previous example in several ways. There are multiple major peaks, distributed throughout the run. There is a long stretch, from about 12 to 21 minutes, where there are many features that may be peaks, or may be artifacts. Close examination of the peak between 24 and 25 minutes will reveal that it is, in fact, at least three separate peaks that are not fully resolved. The major peaks are also taller and wider than the peaks in the previous example.

After examining this file, zoom out (`Ctrl-Z`), and analyze with the default Width and Threshold values. A large number of peaks will be identified. The big peak at the end of the run is not included. The analysis appears to end at around 28 minutes. Zoom in on the area around 27 to 28 minutes. Many peaks which appear to be quite small have been identified and the last peak is number 150. There is a limit of 150 peaks built into the software, which accounts for the end of the data being ignored.


Not all 150 components are of interest. Some parameters need to be changed to reduce the number that are identified. In this part of the exercise, graphical programming will be performed to build the integration table.

Width

Zoom out to full size, then zoom in on the first major peak (around 9.5 minutes). This peak is nearly 0.5 minutes wide at the base, and it is significantly wider than the surrounding smaller peaks. Right click on the electropherogram and select **Graphical Programming**. A new menu of integration parameters will appear. Select Width. A message will appear in the Status bar at the bottom of the instrument window which states "Click on the start of the Peak". Place the mouse pointer at a point near the left edge of the major peak along the baseline and click. You will be prompted to click on the End of the peak. After the Start and End of the peak have been entered, the

following dialog box will appear. Your values may differ slightly from those shown, depending on exactly where you clicked.

Figure 50 Determining Width of Peak



There are several options available in this dialog. Start and Stop times refer to the points where you clicked. For the Width parameter, the value is the difference in these two time points (the width of the peak). There is an option to add this parameter to the Integration Events table or to the Manual Integration Fixes table. These two choices differ in one critical way. Items added into the Integration Events table become part of the method; they will affect any data file that is integrated with that method. The Manual Integration Fixes table becomes part of the data file; items here will affect the integration of only the single data file to which they are attached. In this exercise, we will work only with the Integration Events table. Manual fixes are discussed in the resources in Appendix 1.

The Add to Table option will add the event into the Integration Events table. Analyze Now will add the event to the table, but it will also cause the integration to be redone with the new parameter in place. For now, click Add to Table.

The new Width item will now be visible in the Integration Events table. This item will only be in effect for the interval between the Start Time and Stop Time. We would like this event to be in effect for the entire data set. To have a parameter active over the entire run, set both Start Time and Stop time = 0. Change the Start and Stop times of the new width value to = 0. For purposes of consistency with this manual, change the Value of the new Width event to 0.5 (your value should already be very close to this point.)

Now we have a problem. There are two Width statements with different values, both of which are declared for the entire run. We could delete the old Width statement, but the software offers an easier alternative. At the left of each line is a box with a red check-mark. Click on the red mark on line one to de-select it. This line will be ignored as long as it is unchecked. Your Integration Events Table should look like this one:

Figure 51 Integration Events Table

#	Event	Start Time	Stop Time	Value
1	Width	0.000	0.000	0.2
2	Threshold	0.000	0.000	50
3	Width	0	0	.5
4				

Threshold


Zoom out and click Analyze. You will see a small change in the integration. Obviously, changing the Width alone is not sufficient to integrate this data set. We are still finding too many of the very small peaks that are of no interest. Right click on the electropherogram, select **Graphical Programming**, and then **Threshold**. Following the prompts in the status bar, click near the start of the electropherogram (time 0.0), and at around 8.0 minutes. When the dialog box appears, click Analyze Now. The new threshold value, which should be around 850, will have removed many of the small peaks in the 0-8 minute time range. The new threshold value will also have been added to the Integration Events table. As with the Width parameter, change the start and stop times of the new Threshold event to 0, and deselect the old Threshold parameter. Again for consistency, adjust the new Threshold value to exactly 850. Analyze the sample.

You will observe that the change in threshold has had dramatic effects. Almost all of the small peaks have dropped out of the integration results. The peaks near the end of the electropherogram have been found. But we now have a problem with the cluster of peaks between 24 and 26 minutes. Zoom into this area, and observe that the baseline starts about halfway up the start of the peak. This demonstrates the fact that global parameters do not always work for every peak in the data set. We need to add new lines that affect only the region of the clustered peak.

Integrating a Peak Cluster

Recall that with the original Width and Threshold values (0.2 and 50, respectively) this cluster was integrated nicely. Perhaps one or both of those values, applied to the region around the cluster, will achieve the desired results. While still zoomed in, select Width from the Graphical Programming menu. As prompted, click at about 24 and 25 minutes. In this case, we are using the graphical function here to define Start and Stop times rather than to define Width. In the dialog box that appears, change the Value to 0.2. Click Analyze Now. The cluster should appear similar to the one shown here:

Figure 52 Well integrated peak cluster



This cluster is well integrated. Even the shoulder on the leading edge of the cluster has been detected. Use `Ctrl-Z` to zoom out. There are now two Width values in the Integration Events Table. One is global, the other (the new one) is only in effect over the time range specified in the event.

Integration Off

This integration is nearly done, but there is still a problem at the end of the electropherogram. At the end of the run is a baseline shift that is being recognized as a peak. Zoom into the area from about 34 minutes to the end of the run. We can see that there is really nothing of interest beyond the peak at 35.2 minutes.

To turn off integration, select **Integration Off** from the **Graphical Programming** menu (while still zoomed in). Click at about 35.8 minutes and at the end of the electropherogram, then click **Analyze Now**. The peaks in the specified area will disappear from the integration results, and the **Integration Off** item will be added to the table. Like the **Width** value we used to resolve the peak cluster, the **Integration Off** function is only active over a specific period of time. All of the integration functions except for **Shoulder Sensitivity** can be used globally or locally. For example, you could use another **Integration Off** event from **Start Time=10** to **Stop Time=20** to eliminate the remaining small peaks. For now, the integration appears to be acceptable. The **Integration Events** table should resemble the one below. The **Start** and **Stop** times may differ slightly, depending on where you clicked.

Figure 53 Integration Events table

#	Event	Start Time	Stop Time	Value
1	<input type="checkbox"/> Width	0.000	0.000	0.2
2	<input type="checkbox"/> Threshold	0.000	0.000	50
3	<input checked="" type="checkbox"/> Width	0.000	0.000	.5
4	<input checked="" type="checkbox"/> Threshold	0.000	0.000	850
5	<input checked="" type="checkbox"/> Width	24.322	24.892	0.2
6	<input checked="" type="checkbox"/> Integration Off	35.756	36.992	0
7	<input checked="" type="checkbox"/>			

Delete lines one and two, as they are unused. Leaving them in and unchecked does no harm, however.

Save the method, or all previous efforts will be lost. Select **File | Method | Save As**. Save the method as **Integration Exercise.met** (You may see a warning that the **Separation Time** of the method is 0.0. Select **Yes** to continue with the save. This message appears because, in this exercise, nothing was added to the **Instrument Setup** window of this method, and this omission has been detected by the **Auto-Validation** function of the software. In actual practice, **Integration Events** would be added to the method used to acquire the data.)

Integration Results and Reports

The results of your integration efforts go beyond some labels on the electropherogram. These results may be used to generate reports and to identify and quantitate unknowns. These topics are covered in later sections of this manual. For now, use a pre-defined report to capture the current results. From the menu bar select **Reports | View | Area%**. A report will appear on the screen showing the **Integration** results in tabular form. If a hard-copy of this report is desired (and if a printer connected to your system), right click on the report and select **Print**.

Section 8-Calibration

Introduction

The previous exercise demonstrated how the 32 Karat Software can be used to determine peak parameters such as migration time, area, and height. This information is of limited utility unless it can be translated into terms that are descriptive of the components in the sample. There are two basic questions that may be asked: “what is this component?” and “how much is present in the sample?”

Qualitative Analysis

The first question, “what is this component?” refers to QUALITATIVE analysis. It can be answered in two different ways. First, if a known substance is used as a standard, the migration time or mobility of the resulting peak (possibly combined with other information such as a PDA absorbance spectrum), can be taken as evidence that a peak in an unknown sample represents that known substance. Second, if the migration time or mobility varies in some linear way with a molecular property, such as molecular weight or number of base pairs, a standard curve can be created that can be used to determine the defined quality (such as molecular weight) based on a series of standards. In this method, the unknown need not be present in the standard mixture. Qualitative analysis is covered in the following section of this manual.

Quantitative Analysis

This exercise will deal with the second question, “how much is present in this sample,” which is a question for QUANTITATIVE analysis. In quantitative analysis, we run a series of runs with differing concentrations of a known substance or substances. For each concentration of each component we determine a detector response. By generating a graph of detector response versus concentration, it is possible to determine the concentration of the components in an unknown. The 32 Karat Software simplifies this task. It contains tools for collecting the data from calibration runs, generating appropriate curve fits, and using the curve fit data to analyze unknowns. This exercise will also introduce you to using a sequence for post-run data reprocessing. We will use an External Standard calibration; in this method, data from standards and unknowns are acquired in separate runs.

Creating a Calibration

Included in the Data Samples directory are five files named: Level 1.dat, Level 2.dat, Level 3.dat, Level 4.dat, and Level 5.dat. These data files are the results of runs of the same two components at five different concentrations as indicated in the table below.

Filename	Component (units/mL)	
	<i>Alpha</i>	<i>Beta</i>
Level 1.dat	1.0	4.0
Level 2.dat	2.0	4.5
Level 3.dat	3.0	5.0
Level 4.dat	4.0	5.5
Level 5.dat	5.0	6.0

The first step in the analysis is to identify the expected migration times of the two peaks. To do this, open an instrument window (continue to use the offline UV instrument from the previous exercise). Select **File | Data | Open** and navigate to the Data Samples directory. Open your copy of the file Level 3.dat. Two channels of data are present: Channel A (UV 214nm) and Channel C (current). You may minimize the Channel C window, as it will not be used in this exercise. Channel A is the electropherogram from the analysis of the middle concentration of the calibrators. The smaller of the two peaks between 3.5 and 4.5 minutes is Alpha, and the larger of the two is Beta (previously determined in runs of the individual components).

Developing a Calibration Method

We will build a new method to calibrate our data set. Select **File | Method | Open**. Navigate to your copy of the Data Samples directory and open the method `Calibrate.met`. Select Analyze to integrate the data set. (The integration parameters have already been optimized for you.) The two large peaks will be well integrated with the default parameters.

Right click on the electropherogram and select **Graphical Programming**, then select **Define Peaks**. Click in the electropherogram to bracket the large peaks (approximately 3.7 and 4.2 minutes). The following dialog will appear (your time values may differ slightly from those shown):

Figure 54 Define Peaks dialog

The screenshot shows the 'Define Peaks' dialog box with the following settings:

- Define peaks in range:**
 - Start time: 3.68 Minutes
 - Stop time: 4.21 Minutes
- Migration time window:**
 - Relative: ± 2.5 %
 - Absolute: ± 0.092 Minutes
- Units:**
 - Units: %
 - Quantitate peaks on: Area
 - Minimum peak area: 0
- Actions:**
 - Add all peaks to table
 - Replace existing peaks in table

The Define Peaks function allows you to place any peaks in the specified window into the method's Peaks/Groups table. Only peaks that have been identified during Analysis can be entered into a method. Start and Stop times define the part of the electropherogram that contains the peaks we wish to add. Migration time window allows us to specify how much variability in migration time is allowed before a peak is no longer considered to be the component of interest. This can be defined as either a percentage of the defined time, or as an absolute time interval. We can add all the peaks in the window to an existing table (if any). If we chose to Replace existing peaks, any currently defined peaks that fall within the newly defined time windows will be deleted. For this exercise, accept the defaults.

From the menu bar select **Method | Peaks / Groups**. A table will open listing the two peaks just added. As these peaks have not yet been named, the system will create a name based on the migration time. Click in the Name column and change the name in line 1 to Alpha and the name in line 2 to Beta. Use the scroll bar at the bottom of the window to review the column headings in this table. We will only use some of these columns in this exercise. See the resources in Appendix 1 for more information.

The items that can be changed here include:

- The defined migration time and allowed window for a component
- The data channel to be used (for multichannel detector types)
- The type of curve fit to be used for the calibration data set
- Up to 10 calibration levels


We must input the concentration values from the table above. Scroll right until the columns headed Level 1 through Level 5 are visible. Input the “known” values for Alpha and Beta. Level 1 is the lowest value, level 5 is the highest. For now, we will accept the defaults for the other columns. When all the values have been entered, the table should look like Figure 55.

Figure 55 Peak/Group Tables

#	Name	ID	Level 1	Level 2	Level 3	Level 4	Level 5
1	Alpha	1	1.0	2.0	3.0	4.0	5.0
2	Beta	2	4.0	4.5	5.0	5.5	6.5
3							


Save your method. Right click on the electropherogram and select Annotations. Select parameters so that only “Name” and “ESTD concentration” are in the Show box. Click *OK*. Analyze the sample again. The name of the peak will appear on the electropherogram. ESTD concentration will be 0.000, because we have not yet completed the calibration process. At this point we have only used one of the calibration samples to define our peaks. We will use all five calibrators to generate the standard curves for Alpha and Beta.

Generating the Calibration Curve

To generate the curves, we will reprocess the data through a Sequence. In an earlier exercise, we used a Sequence to perform multiple runs. Sequences can also be used to do batch reprocessing of data. Select **File | Sequence | New** to open the sequence wizard. Click on the file icon  and select `Calibrate.met` (the file we have been editing). For Data File Type, click to select From existing data files. Click *Next*. The next dialog allows us to select multiple data files. Click on the file icon. This will open a dialog that allows the opening of multiple files at one time. Navigate to the Data Samples folder. Locate the file `Level 1.dat`. Double click on the file name. When you double click, the file name is added to the list at the bottom of the dialog. Repeat for Level 2, Level 3, Level 4, and Level 5, so that all five files appear in the list. Click *Open* when done. Verify that the five files are listed in the window, then click *Finish*. This will open a Sequence table with the method and the five data files listed. Scroll

right to the Filename column. The five files should be listed in order. Scroll back left to the Level column. Change the level in line 1 from 0 to 1. Use the down arrow to move to the same column in line 2, and change the 0 to 2. The Level entered in this column of the sequence table corresponds to the Level column in the Peak \ Group table. Continue down the column changing the remaining zeros to 3, 4, and 5, respectively. As you change the level to a non-zero value, the sample type automatically changes to Calibration. Click in the Run Type column of line 1. The dialog shown in Figure 56 will open. The items here are the available run types. Click on Clear All Calibration and click OK. The run type for line one will change to “CAL CCA”. Selecting this type will cause any previous calibrations of this data set to be cleared before the new calibration is performed.

Figure 56 Sample Run Types dialog



Save the sequence with the name `Calibrate.seq`. The right and left sides of the window should look like Figure 57.

Figure 57 Sequence Table

The figure displays two overlapping screenshots of the 'Sequence: Calibrate.seq' window. The top screenshot shows a table with the following data:


Run #	Method	Filename	Sample Amt.	ISTD Amt.	Mul
1	Calibrate.met	level 1.dat	1	1	
2	Calibrate.met	level 2.dat	1	1	
3	Calibrate.met	level 3.dat	1	1	
4	Calibrate.met	level 4.dat	1	1	
5	Calibrate.met	level 5.dat	1	1	
6					

The bottom screenshot shows a table with the following data:

Run #	Status	Run Type	Level	Conc Override	Reps
1		Calibration	1		1
2		Calibration	2		1
3		Calibration	3		1
4		Calibration	4		1
5		Calibration	5		1
6					

We are now ready to Process the sequence and generate the calibration. Select Process from the Sequence menu. The dialog in Figure 58 will appear. Run range must be All. Processing mode must be Reintegrate. Click Start. The processing will require a few seconds to complete. As each data file is analyzed, the status column in the Sequence table will change to “Complete”.

Figure 58 Process Sequence dialog



From the menu bar select **Method | Review Calibration**. A window will open showing the results of the calibration for the first component (Alpha). The data looks reasonable, but the line fit could be better. Select Beta in the upper right corner of the window. The curve fit here is certainly not acceptable. The lower right corner of the window indicates that this fit has been made point-to-point. Perhaps another type fit would be better? Re-open the Peaks / Groups table. Scroll left to the column headed “Fit Type”. For line one (Alpha) select Quadratic. For line 2 (Beta) select Linear. In line one, click the box in the Zero column to force the line through the origin. Do not click Zero in line 2. Save the method and look at the calibration results again. This time the curve fit should be better for both Alpha and Beta. In practice, you would select and vary the curve fit parameters to determine the most appropriate fit for your data set.

Saving the method will save the calibration. It is not necessary to reprocess the sequence.

Analyzing Unknowns

Creating a standard curve is only useful if it can be used to determine the amount of material in an unknown. Unknowns can be added to the sequence table after calibration lines. When processed, the unknown peaks will be compared to the standard curve and the values calculated. Open the sequence window. In line 6, click in the Method column and select `Calibrate.met` (the method we have been using). Under Filename, select `Unknown1.dat`. Save the sequence. Instead of reprocessing the entire sequence, click on the Range button in the Process Sequence dialog (Figure 58) and enter 6 as the line to be run. Click Start. To see the results of this analysis, you can display the electropherogram for `Unknown1.dat`, or you can view a report. To use a default report, select **Reports | View | External Standard** from the menu bar. To print this report to your system's default printer, right click in the report window and select Print.

Creating custom reports is covered later in this manual.

This exercise has covered only the basics of quantitative analysis. For more detail and more options, refer to the resources in Appendix 1.


Section 9-Qualitative Analysis

Data Analysis and Reporting

In some modes of CE separation, the order of migration is determined by a molecular quality. The most common applications that utilize this approach are separation of proteins by molecular weight and the separation of nucleic acids by number of base pairs. The 32 Karat Software allows you to construct a standard curve based on parameters such as these. In the case of DNA, for example, a series of fragments of known length are separated in a gel-filled capillary. A plot of migration time versus chain length is constructed, and this plot is used to determine the length of an unknown fragment.

In this exercise we will create a standard curve based on the base number of single stranded DNA. A sample of ssDNA is injected into a capillary containing a polymer network. Smaller molecules can negotiate through the polymer strands more easily than can larger fragments, resulting in a separation based on DNA chain length. The migration time of peaks of known lengths of ssDNA is used to create a standard curve that can be used to estimate the length of an unknown fragment,


Figure 59 Electropherogram of Single-stranded DNA



Open an offline instrument configured for a UV detector. From the Data Samples directory, open the data file `ssDNA_40-60.dat` and the method `Qualitate.met`. The data file represents the separation of a mixture of poly-A oligonucleotides ranging from 40 to 60 bases in length. The method contains an Integration Events table that will successfully integrate these peaks, and a Peaks \ Groups table that will name the peaks with the number of bases they represent. Analyze the data. Right click on the electropherogram and select **Annotations**. Add Name to the list of annotations to display, then click OK.

Next, we must enter the base numbers (Y) and the migration times (X) into the method. From the method menu select "Qualitative analysis" to open the dialog shown below.

Figure 60 Qualitative Analysis dialog



We can select the type of data we want for the X-axis. In this exercise, we are using migration time, but we might also use another parameter such as Mobility. We can select either a linear or a logarithmic scale for this axis. Minimum and Maximum refer to the X-axis limits over which we will allow the interpretation of qualitative data. This limits how far beyond the ends of the data set we are willing to extrapolate to find the quality of an unknown.

For the Y-axis, we can use any name we chose. Because our quality is the number of bases comprising the oligonucleotide, enter Bases. A Linear Y-axis scale will be appropriate for our data.

There is a two column spreadsheet in this dialog box (Figure 60). The Left column contains data for the Y-axis. The heading of this column will contain the text entered for the Y-axis label, in this case “Bases”. The Right column contains data for the X-axis. The heading is that selected for the X-axis. In this case we are using Migration Time. You must manually enter the X and Y values into the spreadsheet. One way to obtain this data in tabular form is to open an Area % report by selecting **Reports | View | Area %** from the menu bar. This may be printed if you have a printer connected to the system, or the report can be displayed as a window alongside the Qualitative Analysis window. For each peak, enter the base number (from 40-60) and the corresponding migration time. As you enter the data, it will be graphically displayed in the lower right corner of the dialog box.

When all the data have been entered, select the Fit Type to determine the best curve fit to the data set. In this example, the highest “Goodness of Fit” values is obtained using a Cubic fit. When done, do a Save As; rename the method `Qualitate 1.met`.

Analyze the data file again. Right click on the electropherogram and select **Annotations**. Add “Quality” to the right column (“Name” should already be present). Click on the word Quality to highlight it, then set the number of decimals to 2. Click OK. The electropherogram will now display the assigned base number (as Name) and the base number calculated using the standard curve (as Quality). The base number is an integer; the small discrepancy between the assigned value and the calculated value is due to the fact that our goodness of fit was not perfect.

This method can now be used to estimate the size of an oligonucleotide of unknown length. To do this, the unknown would be analyzed by the same method used to analyze these standards, and the data analyzed with this method.

Section 10-Using Mobility

Mobility is a parameter that quantifies how a charged particle migrates in an electrical field. The equations for calculating mobility are defined in Appendix 2. Stated briefly, a component with a higher mobility will move more rapidly through the separation medium than will a component with a lower mobility. Because particles may be attracted to either the cathode or the anode, mobility has a vector component — it may be positive or negative. Mobility toward the cathode (negatively charged electrode) is defined as positive, and mobility toward the anode (positively charged electrode) is defined as negative.

Mobility is not a constant like molecular weight, as it can change depending on the separation conditions employed. As an example, consider the amino acid glycine. At alkaline pH, this molecule has a net negative charge ($\text{NH}_2\text{-CH}_2\text{-COO}^-$) and a negative mobility. At acidic pH, it has a net positive charge ($\text{NH}_3^+\text{-CH}_2\text{-COOH}$) and a positive mobility. Near pH=7, it has no net charge ($\text{NH}_3^+\text{-CH}_2\text{-COO}^-$) and a mobility of zero.

In practice, mobility is defined for a given molecule under a given set of conditions. Variations in separation conditions that affect all species equally (such as variations in electroosmotic flow or the voltage delivered by the power supply), and variations that have similar effects on closely related groups of compounds (such as the effect of small pH changes on a series of basic drugs) can be accounted for by including in the analytical run a standard of defined mobility.

There are two different mobility values to consider: electrophoretic mobility (μ), which represents the movement of the molecules in the electrical field, and apparent mobility (μ_{app}), which is the algebraic sum of the electrophoretic mobility and any other forces which drive the sample through the capillary. The most common additional force is electroosmotic flow (EOF, sometimes referred to as μ_{eof}), although other factors such as pressure-assisted migration could also be considered here. In this discussion, we will limit apparent mobility to the sum of electrophoretic mobility and EOF:

$$\mu = \mu_{\text{app}} - \mu_{\text{eof}}$$

The 32 Karat Software has the capability to use mobility as an alternate to migration time for peak identification. Used in this way, mobility increases the robustness of a capillary electrophoretic method by correcting for small changes that may be present from run to run. The software is also capable of displaying electropherograms with the X-axis scaled in terms of mobility rather than time. In this exercise, you will review these features of the software.

Mobility Markers


Every run that will incorporate mobility calculations must include a mobility marker. This is a component for which the mobility has been defined under a given set of run conditions. Reproducing these run conditions in subsequent runs is critical in using mobility as a parameter for peak identification. Buffer factors, such as ionic strength and pH, are especially important. See the resources in Appendix 1 for more information on this topic.

Assigning Mobility Values

Assigning a mobility value to the mobility marker involves a series of steps. The following assumes that the method development work to create an effective separation technique for the samples to be analyzed has already been developed.

1. Create a sample containing the proposed marker and an EOF marker. The sample matrix should be a close approximation of the unknown samples to be analyzed. The EOF marker is a small, uncharged molecule that can be detected at low concentration. It is used to calculate the contribution of electroosmotic flow (μ_{eof}) to the apparent peak mobility. If running under conditions where EOF is minimal (such as a neutral coated capillary) the EOF marker can be omitted, as $\mu_{\text{eof}} = 0$. Standards for other peaks expected in the unknowns should be included in this run.

- Run the separation, and integrate the results. Adjust all integration parameters to their final form. Add the mobility marker peak and the EOF marker peak (if used) to the Peaks / Groups table. The method must contain the correct dimensions for the capillary length. These are entered in the Capillary/Performance tab of the Advance Methods Options window (accessible from the Method menu). The capillary length and length to the detector are required.




- Annotate the electropherogram with Apparent Mobility and Mobility. Analyze the data again. The Apparent Mobility and Mobility values will be displayed. At this point, the mobility values should all be zero.
- Calculate the electrophoretic mobility of the mobility marker according to:

$$\mu = \mu_{\text{app}} - \mu_{\text{eof}}$$

where μ_{app} is the apparent mobility of the reference peak and μ_{eof} is the apparent mobility of the EOF marker.

- Scroll right in the Peaks / Groups table to the Mobility column. In the row containing the mobility marker, enter the value of μ determined in step 4. Click the μ box in the Mobility Marker column of that same row. Save the method.

6. Analyze the sample again, and examine the electropherogram. The mobility value of the marker should now be non-zero, as should the mobility values for all other peaks in the electropherogram except the EOF marker. The EOF marker should have a value very close to zero (variation from zero at this stage is the result of rounding errors).
7. If the results are acceptable, the method can now be edited to allow the generation of an electropherogram with mobility as the X-axis in place of time. This option is selected in the Initial Conditions tab of the Instrument Setup window:



Mobility Plot

Selecting Mobility (or Apparent Mobility) will cause the creation of an additional plot after data analysis. This plot will show the electropherogram re-scaled so that the X-axis is displayed in terms of Mobility (or Apparent Mobility) rather than time. The “Plot trace after voltage ramp” option should be selected unless peaks are detected during the voltage ramp-up phase of the separation. Deselecting this option may result in the voltage ramp interval being over-emphasized in the resulting mobility plot.



NOTE

Mobility channels cannot be generated for data files that were created with this option turned off. It must be ON during data acquisition. When using this option, there must be at least one Mobility Marker identified in the Peaks/Groups table of the method.

After selecting the appropriate mobility channel options, save the method. Run the sample again. At the end of the run the mobility trace will be generated.


Demonstration: Using the Mobility Features



NOTE

The following is not intended to be a hands-on exercise, but rather a step-by-step example of the process described above. If you wish, you can attempt to reproduce these results on your instrument. This example utilizes a UV detector and Beckman Coulter Test Mix B (P/N 501333).

We have previously determined the optimum conditions for the separation of a mixture of two components, ALPHA and BETA (see the exercises earlier in this manual). ALPHA has been chosen as our mobility marker. Our method has considerable EOF, so we will need to use an EOF marker. We have chosen to use Ethanol; it is readily soluble in our sample, creates a detectable peak, and does not interfere with our analysis (all previously determined -- in practice, this must be shown by experiment). A final concentration of 5% Ethanol has been chosen. The sample is prepared and run. After Step 3 (above) the resulting electropherogram looks like this:



Our ethanol marker gives a small, negative peak, but it is properly detected by the software. The values obtained are listed in the table below:

Table 7 Mobility values, first pass


Peak Name	Apparent Mobility	Mobility
EOF	0.00060227	0.00000000
ALPHA (marker)	0.00034763	0.00000000
BETA	0.00032093	0.00000000

As expected, the faster moving EOF marker has the higher apparent mobility. The units of mobility are $\text{cm}^2\text{V}^{-1}\text{sec}^{-1}$


We can calculate the Electrophoretic Mobility of ALPHA according to

$$\mu = \mu_{\text{app}} - \mu_{\text{eof}}$$

substituting yields $\mu = 0.00034763 - 0.00060227 = -0.00025464$, the mobility of ALPHA under these separation conditions. The mobility value is a negative number, indicating that ALPHA's electrophoretic migration is toward the anode, and is opposite the direction of the electroosmotic flow.

 NOTE	<p>In practice, the above run should be repeated several times to verify that the system is giving consistent results before defining the mobility value for the mobility marker.</p>
---	---

We now re-open the Peaks / Groups table and insert the value of the mobility marker as described in Step 5. The table now looks like this.



#	Name	ID	Mobility Marker	Mobility	Mobility Window	
1	EOF marker	1	<input type="checkbox"/>	0.00000000	0.00000000	PDA
2	ALPHA	2	<input checked="" type="checkbox"/>	-0.00025464	0.00000000	PDA
3	BETA	3	<input type="checkbox"/>	0.00000000	0.00000000	PDA
4			<input type="checkbox"/>			


We save the method and analyze the data file again. The results are now:

Table 8 Mobility values, second pass

Peak Name	Apparent Mobility	Mobility
EOF	0.00060227	-0.00000000
ALPHA (marker)	0.00034763	-0.00025464
BETA	0.00032093	-0.00028134

The electrophoretic mobility of the EOF marker is close to zero, as expected. The mobility of ALPHA is the value we defined it to be. The mobility of BETA has now been calculated relative to the mobility of ALPHA.

We can now set up to collect data for a mobility plot, if we wish. Open the Instrument Setup window and click on the Initial Conditions tab. We select Mobility (see the figure above). After saving the method and running the sample again, we have two electropherograms. One is the UV electropherogram, which looks very similar to the run we used to measure the mobility of the marker. The second appears VERY different, as shown here.




The order of the peaks has been reversed, because the latest eluting peaks have the most negative mobility. The relative spacing between the peaks has also changed.

Peak Identification by Mobility

In the exercise on mobility, we used migration time and a migration time window to establish the identity of peaks in a data file. In the same fashion, we can use mobility and a mobility window to identify peaks, based on the mobilities of standards. The choice of migration time or mobility applies to the entire electropherogram; it is not possible to use time for some peaks and mobility for others.

The option to use mobility for peak identification must be selected on the Options tab of the Method Properties window (under the Method menu). Select “Mobility based” and click OK. Next, peaks are added to the Peaks / Groups table (if not already present). This is most easily performed by using the “Add Peaks” function from the Graphical Programming menu as described earlier in the section on Integration. Open the Peaks \ Groups table and scroll to the Mobility columns. The assignment of the mobility of a mobility marker has already been described. To identify other peaks by mobility, it is necessary to manually input their mobility value and window. For the peak BETA in our example, the mobility value was determined to be -0.00028134, so this value is entered into the table for the line defining BETA. The default Mobility Window is zero, which means that the peak would only be detected if it had exactly the same mobility value as that defined. Some variation can be expected. The appropriate value will be determined by trial and error. A window of 5-10% is a good starting point. For BETA, a 5% window would be 0.00001407; the window is always an unsigned number (not positive or negative). At this point our Peaks / Groups table looks like this:



#	Name	ID	Mobility Marker	Mobility	Mobility Window	
1	EOF marker	1	<input type="checkbox"/>	0.00000000	0.00000000	PDA
2	ALPHA	2	<input checked="" type="checkbox"/>	-0.00025464	0.00000000	PDA
3	BETA	3	<input type="checkbox"/>	-0.00028134	0.00001407	PDA
4			<input type="checkbox"/>			

When we click on Analyze, the peak identification will be by mobility. To verify this, change the Mobility entry for BETA to some other value like -0.00025227 and Analyze again. Because the mobility of the peak is now outside the window, it will be integrated, but not associated with the ID “BETA”.

In the section on Reports we will create a report template that will contain mobility data.


Section 11-Creating Reports

Introduction

The ultimate purpose of operating a CE or any other instrument system is to answer a question. This answer is only useful if it can be communicated. The 32 Karat Software includes a powerful suite of utilities for creating printed reports. It also allows you to export data and results to other software packages. Full integration into the Windows NT environment allows you to copy and paste electropherograms and other information from 32 Karat Software directly into a wide variety of other software packages, such as word processors. This section of the manual will describe how to use the built-in report generating features of the software. For information on other ways of moving data to other applications, see the resources listed in Appendix 1.

In other sections of this manual, you have been introduced to some of the default reporting options. Here you will learn how to customize a default report to include the information you want to present. We will continue to use a UV instrument as a prototype. LIF reports are essentially identical to UV reports. PDA data can be reported in the same formats as UV data; additionally, other PDA specific report types are available such as PDA Spectrum Library reports.

The most simple approach to reporting is to use pre-defined reports. These were used in earlier sections of this manual. From the menu bar, select **Reports | View** to display the menu shown below which lists the default set of pre-defined reports.



In this exercise we will use the data file and method from the exercise on Qualitative Analysis. A review of the reports menu will show that there is no default report for qualitative data, so one will have to be created. We will create a Method Custom Report, (which will become part of the method). We will also save the report as a standard template that will be added to the reports list.



Creating a Method Custom Report

Open your Virtual UV instrument from the main 32 Karat window. From your working folder, open the file `Data Sample 3.dat` and the method `Quality.met`. Use the copies of the data and the method that were used during the exercise on Qualitative Analysis; we will be creating a report based on the analysis results generated during that activity.

From the **Method** menu, select **Custom Report**. A blank window will open. (If a method custom report has previously been developed with the method you are using, it will automatically open at this point.)

Modifying a Report Template


Rather than developing a report from nothing, we will modify an existing report template to simplify this task. A report template is a pre-defined set of instructions for the creation of a report. A number of templates are contained in the software. You can add templates by creating a new one or editing an existing one and saving it under a new name. We will use the `Area%` template as the basis for our Qualitative Analysis report. From the file menu select **Report Template | Open** and select `Area%.srp`. The template will open with the results from the open data file displayed. The indications that we are looking at a template and not at a report are that the header information is highlighted, the table has a broken border, and the ruler is visible.

As a first step, we will change the title of the report. The report editor works very much like a word processor. Using the mouse, highlight the heading "Area % Report." Type `Qualitative Analysis Report`. Highlight the new title, and select **Bold**  and **Italic**  from the toolbar.

The information below the title that is highlighted in gray represent fields. Fields allow you to insert blocks of information without tedious code writing. Fields and text may be used together; in this template, each field is preceded by a text descriptor. These fields may be easily repositioned. Click after the "Acquired" field, and press the Delete key. This will move the Print Time field to the same line as the Acquired field. Press the Tab key twice to place space between the two fields. Click in the blank line below these two fields and type `User ID:`, then press the Tab key. Right click and select **Insert Field**, then **User Name**. A new field will appear with the current User Name (the default name "System" will appear if no user name has been defined). Use the keyboard just as you would a word processor to align the field and the text in the desired manner.

Below the header information is a graph displaying the electropherogram. By default, the graph is scaled to show the full range of data. In the case of Data Sample 3, the area of interest is only at the right end of the electropherogram. As the next step, we

will adjust the electropherogram display so that the peak region is more easily seen. Double click on the graph to display the Graph Properties box.




The Trace Setup tab is where we select the data to be displayed. We can select the Current Data (the file currently open in the Instrument window) or a data file (or files) stored on the computer. For example, we may want to display the results from analysis of a standard along with results from an unknown. We can display multiple channels from a common data file by selecting the same Data Source, but selecting different channels in the Trace column. For example, the UV and electrical Current traces could be displayed simultaneously. The X and Y offset and scale features can be used to position the traces so that the pertinent details are visible. These features are described in the resources listed in Appendix 1. Units is a label which can be adjusted to match the data displayed. For an LIF detector, this could be changed to “Fluorescence” by clicking in the Units box and entering the desired text.

Click on the Annotations button. Because this will be a Qualitative report, remove any existing annotations from the right side panel. Add Quality, and set the number of decimal places to zero. Click OK to return to the main dialog.

Now we can change the X-axis so that the peak area fills the graph window. Click on the Axis Setup tab. With X-axis selected, click on Use this range. Enter a minimum of 24 and a maximum of 30 to bracket the 24-30 minute time interval. When done click


OK. The part of the electropherogram containing the 21 peaks should fill the graph, and each peak should be properly labeled with the Quality value.

The graph can now be resized by dragging the handles on the right and bottom edges of the image. Moving the graph requires that it be treated as a text object. To center it on the page, single click the graph to select it (a narrow dotted border will appear) and use the  icon from the toolbar. Space can be added above and below the graph by inserting blank lines with the Enter key.

Additional charts can be inserted by right clicking on an empty line and selecting **Insert Graph**. You might, for example, want a second graph showing the full time range of the data set. The new graph is then formatted as described above.

Below the graph is a table with a variety of peak parameters. Some of these we will keep, others we will change for our new report. Right click on the column header "Area%" and choose **Change Parameter**. Select **Quality** as the new parameter, and change Decimals to one. In the Column Header box type Bases. Click OK. We now have a column that contains the number of bases that the software has calculated for each peak.

Right click in the table and select **Report Properties**. We are not interested in Height or Height %, so remove them from the right-side box using the red arrow. Click OK.

Our table is now too wide for the remaining three columns. By clicking on the vertical dividers between columns, we can drag the column widths to whatever we desire. At any point during editing, we can click on the Print Preview icon  to see how the report will look on our default printer.

Save the method. The new report is now part of the Quality method, and will be available whenever the method is used by selecting Method Custom Report.

Saving a Report as a Template

In order to use the report in another method, we can save the it as a template. From the File menu, select **Report Template | Save As**. Type `Quality` as the filename and click. By default, the template will be saved with the extension `*.rep`. Select **File | New Method**. Open Method Custom Report by selecting it from the Method menu (if not already open). The report will be blank, because a new method has no defined Custom Report. From the File menu, select **Report Template | Open**, and open the `Quality.rep` template that you created. When the new method is saved, it will contain the `Quality.rep` as the default Method Custom Report. The Custom Report saved in a method can be edited and changed at any time.

We can also save this report as one of the default reports in the Reports menu. Again select **Report Templates | Save As** from the File menu. This time save the template with the name `Quality.srp`. Click Save. By saving with the `*.srp` extension, you have created a new default template. The next time an instrument window is opened, the new report type will be available from the **Reports** menu. This is a handy way to handle reports that are used frequently.

Appendix 1

Additional Resources

32 Karat Online Help

The Online Help is your primary software reference for the 32 Karat Software. It describes all of the features used in this manual. In addition, it describes all of the advanced features which are not covered here. To access the Help, select **Help | Contents** from the instrument window menu bar. Context-sensitive help is available by selecting the Help button located in dialog boxes, or by pressing the **F1** key.

Installation and Maintenance Manual

The first place to look for answers to questions on hardware and maintenance is the Installation and Maintenance Manual. This document includes detailed descriptions of routine operations, such as rebuilding a capillary cartridge, and of non-routine operations such as minor repairs. A printed copy of this manual is shipped with the instrument. An electronic copy is available on the Manuals CD-ROM shipped with the 32 Karat Software.

Beckman Coulter Service

If you have determined that a problem exists with the instrument that can not be fixed by following the procedures in the Installation and Maintenance Manual, contact your local Beckman Coulter, Inc. Service Representative. There may be a charge for this service.

Help by e-mail

A special mailbox has been created for users of 32 Karat Software. Send your questions and comments to 32karat@beckmancoulter.com. A team of specialists will review your message and contact you. Messages received outside of regular business hours (Monday - Friday, 8:00 am to 5:00 pm Pacific Time) will be reviewed the next business day.

Other references

There are many publications dealing with applications of capillary electrophoresis, methodology, routine operations, and related techniques. Your local research library is the best place to begin. The following list is not intended to be all inclusive, and inclusion here does not imply endorsement by Beckman Coulter, Inc.

Beckman Coulter web site. <http://www.beckmancoulter.com/> - *The online source for the latest information on CE from Beckman Coulter.*

P/ACE Setter Newsletter - *A quarterly newsletter for scientists employing capillary electrophoresis. P/ACE Setter features tips and techniques and product updates. Each issue features a specific application for CE. Free subscriptions are available through the Beckman Coulter web site, or contact your local Beckman Coulter representative.*

Your Local Beckman Coulter Representative - *The best source for product information. Contact your representative for information on the latest products, or to place an order.*

Appendix 2

CE Calculations

Processing of electrophoretic data is similar to the processing of chromatographic data, but there are some differences. Some “CE Specific” calculations are built into the 32 Karat Software. This appendix describes those, along with some other useful calculations.

Ohm's Law

$$V = IR$$

The most important relationship to know in capillary electrophoresis is Ohm's law. This simple relationship defines the interaction between Voltage (V), Current (I), and Resistance (R). For example, if Voltage is held constant, change in Current must be due to a change in Resistance. Information of this type is invaluable in developing and troubleshooting CE methods.

Corrected Migration Time

$$t_{corr} = t \left(\frac{t'_{ref}}{t_{ref}} \right)$$

t = migration time of peak of interest

t_{ref} = actual migration time of reference peak

t'_{ref} = expected migration time of reference peak (from peak id table)

Corrected Area

$$A_{corr} = vA = \frac{L_d A}{t}$$

v = velocity

L_d = capillary length to detector

t = migration time

A = uncorrected peak area

In CE, early peaks migrate through the detector window faster than do later peaks. This creates a peak area bias, which is eliminated by using corrected peak area.

Mobility

Mobility is calculated as the apparent mobility minus the contribution of electroosmotic flow.

$$\mu = \mu_{app} - \mu_{eof}$$

Once the mobility of a reference analyte has been determined, the mobility of related analytes can be calculated.

$$\mu = L_d L_t \left(\frac{1}{Vt} - \frac{1}{V_{ref} t_{ref}} \right) + \mu_{ref}$$

V = average applied voltage up to the migration time of the peak of interest

L_d = capillary length to detector

L_t = total capillary length

t_{ref} = migration time of reference peak in the current run

μ_{ref} = defined mobility for the reference peak

V_{ref} = average applied voltage up to migration time of reference peak

t = migration time of the peak of interest

Apparent Mobility is the sum of true mobility and the mobility caused by the electroosmotic flow. Apparent mobility can be calculated directly for any peak. The apparent mobility of a neutral marker (one that has a true mobility of zero) is called μ_{eof} , the electroosmotic mobility.

$$\mu_{app} = \frac{v}{E} = \frac{L_d L_t}{Vt} \text{ (apparent mobility)}$$

L_d = capillary length to detector

L_t = capillary total length

V = average voltage up to migration time of peak

t = migration time of the peak

v = velocity

Mobility calculations depend on an accurate measurement of the average applied voltage up to the point that the component passes the detector. The output from the P/ACE MDQ voltage monitor is used to calculate average applied voltage as follows:

$$V = \frac{\sum_{i=1}^n V_i}{n}$$

n = data point number at peak migration time

V_i = voltage value at data point i

Appendix 3

Understanding PDA Data

Most of the exercises in the manual have used example data from an ultraviolet detector system. LIF data files appear nearly the same as UV data files. From the point of view of data analysis, the Y-axis scale is the primary difference between UV data and LIF data.

PDA data consists of multiple components. These are:

Channel data: each channel is equivalent to a UV trace. Up to three channels can be collected simultaneously.


3D spectral data: the three axes of this data set are time, wavelength, and absorbance.

Multi-Egram data: similar to channel data, this is a trace of absorbance versus time extracted from the 3D spectral data. There are no limits to the number of data sets that can be created.

Channel data and multi-egram data are processed in the same manner as UV data; all the procedures described in this manual apply.

The remainder of this appendix will describe some of the features of PDA data. If you wish to re-create these displays on your computer screen, you will need to open an instrument configured for a PDA detector (offline mode is satisfactory). From a copy of the Data Samples directory, open the data file `PDA Data Sample.dat` and the method `PDA.met`. In the figures in this manual, the background color of the windows have been changed from the default color scheme to facilitate printing.


The multiple components of the PDA data set are available by choosing **View | Gallery** from the menu bar.



Each item in the menu represents a different way of presenting the data on the screen. Only one of these selections may be made at one time.

Mixed View w/ 3D option

This option is shown below. This image illustrates several of the other options as well.



This window shows four different views of the same data set. The sample is the same mixture of Alpha and Beta that has been used throughout this manual. This window is divided into four panes. By clicking and dragging the bars that separate the panes, the relative areas assigned to each pane can be changed. The views available here are:

	Left	Right
Upper	Contour Plot	Spectrum
Lower	Electropherogram	3D Plot

Contour plot

The Contour Plot is an overhead view of a three dimensional data set. The X and Y axes are Time and Wavelength, respectively. Absorbance at any given wavelength and at any given time is indicated by color or by a shade of gray. This plot can be “zoomed” by using the mouse to “rubber box” the desired area. Right click and select **Full Unzoom** to restore the view.

Spectrum


The spectrum detected at a specified point in time. This is selected by moving the horizontal slider in the Contour plot to the desired time point. It represents a slice through the 3D data. This plot can be “zoomed” by using the mouse to “rubber box” the desired area. Right click and select **Full Unzoom** to restore the view.

Electropherogram

The electropherogram at a specified wavelength. This is selected by moving the vertical slider in the Contour plot to the desired wavelength. It represents a slice through the 3D data in a direction that is perpendicular to the spectrum view. This plot can be “zoomed” by using the mouse to “rubber box” the desired area. Right click and select **Full Unzoom** to restore the view.

3D Plot

A three-dimensional representation of the data. The presentation of this window is not tied to the contour view. To adjust the display, right click on the 3D plot and select properties to open this dialog box:



The plot can be displayed as a one-color surface or color can be used to specify absorbance ranges. Colors are user selectable. The image can be rotated around the X and Y axes by clicking on the arrow buttons or by typing in Elevation and Rotation parameters. Zooming is accomplished by setting limits of time, wavelength, and absorbance. The Apply button causes the changes to be made to the 3D plot without closing the properties window. This allows you to test changes without having to close and re-open the dialog.

Mixed View is the same as shown above, except that the 3D plot is not displayed. This option is handy for speeding up the display of very large data files.

3D and Contour Plot shows only these two plots in a two-pane window

The individual **3D**, **Contour**, **Electropherogram**, and **Spectrum** plots display the individual pane as a full window. The limits of the Electropherogram and Spectrum plots are those set on the Contour window, even if the Contour window is not being displayed. The remaining choices in the Gallery menu are not related to these four items.

The **Max Plot** view is an electropherogram that does not represent a particular wavelength. It is built by taking the maximum absorbance across the wavelength scan range at each point in time. This plot is useful for finding peaks that might otherwise go undetected. This data may be analyzed and the results reported.

The **Multi Electropherogram** view is only available if multi-electropherograms are defined in the Instrument Setup component of the method. Like the Max Plot, these traces are synthesized from the 3D data set. To set up these channels, open the Instrument Setup dialog form the Method window and click on the Multi-Egram tab. Three items must be specified for each channel. The channel can be checked On or Off. Channels which are Off will not be displayed and cannot be analyzed. Wavelength refers to the central wavelength of the data to be incorporated into the multi-electropherogram channel. Bandwidth specifies how many adjacent wavelength bands will be combined to form the new data set. These bands are weighted, with bands more distant from the central wavelength contributing less to the final result than do bands closer to the midpoint. Using too narrow a bandwidth will result in a noisy data file. Using too wide a bandwidth will cause a loss of spectral sensitivity in the data set.

These data may be analyzed and the results reported. If many Multi Electropherogram channels have been defined, each will occupy a small area of the display. The frame around each channel can be moved by clicking and dragging with the mouse. Channels can also be temporarily disabled by unchecking them in the Multi Egram tab of the Instrument Setup Dialog.

Ratio opens a three pane window. Two windows will display the electropherograms at specified wavelength. The third displays an electropherogram that is the arithmetic ratio of the two ($Abs\#1/Abs\#2$, on a point by point basis). The parameters for this window are set on the X:Y tab of the PDA Setup dialog (available from the Method menu).

When PDA data is selected, the drop-down box in the tool bar will contain a list of the available data channels. Before data can be integrated or analyzed, the desired channel must be selected from this list. Each data channel will have a unique Peaks \ Groups table and Integration Events table, that will be available from the Method menu only when the data channel has been selected.

More information about using PDA data is available in the resources in **Appendix 1**.

Appendix 4

System Administration

The System Administration functions include adding and configuring instruments, which has already been discussed in this manual. It also includes features for defining users and projects. Projects definitions include: defining default folders in which data and methods will be stored; specifying what instrument(s) can be used for a project; defining which users can work on a project, and what privilege level each user will have. These administrative functions are handled through software “wizards,” which are pre-defined step-by-step processes. It is NOT necessary to implement this function in order to use the system.

System Administration is managed from the main 32 Karat Software Screen. From the Tools menu select Options. Click on the Enterprise tab. Clicking on Enable user login and project management will activate the system administration function. The effects will be seen the next time the program is started.

If the System Administration function is activated, a system administrator must be defined immediately. If no system administrator has been defined, there will be no user with the privilege to disable the system administration function, should that be necessary. Depending on the privilege levels assigned to users, not having a system administrator may make some or all software features inaccessible.

When the Enable option is selected, the User list becomes available. Clicking on Add user opens a dialog in which the user's name and password are entered. A password is not required. If a password is used, it must be recorded elsewhere, as there is no feature in the software for generating a list of user passwords.

When all the names are entered (or at least the name of the administrator), user privileges must be assigned. A user with no privileges cannot use the software! New users do not have any privileges until privileges are assigned.

Click *OK* to close the Options dialog. From the Tools menu select System Administration Wizard. The wizard will only be available if System Administration Mode has been selected (checked). A screen will open with three options for: Users, Instruments, and Projects. If there are secondary administrators, they may be given access only to some of these wizards.

The **Users Wizard** allows the assignment of privileges. Follow the wizard step-by-step. Each user can be given access to specific areas of the software and to specific instruments. Selection will depend on the functions of the individual and the degree of security required in the operation of your instrument systems. For example, regulated quality control laboratories generally implement tighter controls over instrument and data access than do laboratories engaged in basic research.

The **Instrument Wizard** allows the administrator to assign users to specific instruments, but not to assign specific privileges on those instruments.

The **Project Wizard** allows the administrator to create projects. A project is a combination of users, instruments, and folders. Projects are most valuable when instrument systems are used by multiple operators or groups who do not need access to one another's methods and data. Creating a project allows the operations on the instrument to be compartmentalized.

The software is shipped with one project, which is called Default. This is a special project, in that it cannot be protected. It should not be used for methods and data that are not intended to be shared.

For more information on System Administration, see the resources in Appendix 1.